

Оборудование для ретроградной холангиопанкреатографии AXS_tome, easyCure, GPS Wire, twist'n'CATCH, ROCKSTAR, KENA Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48

Россия +7(495)268-04-70

Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81
Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73

Киргизия +996(312)-96-26-47

Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35

Казахстан +7(7172)727-132

Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93



Sphincterotomes

Smooth cannulation performance by Micro-Tip⁺-Design

Precise tip control

Safe guide wire exchange


AXS_tome⁺

When precision matters most

AXS_tome⁺ offers you reliable orientation as well as soft and gentle probing thanks to its improved Micro-Tip⁺ design.


A user-friendly short wire port and optimized X-ray visibility make AXS_tome⁺ the ultimate in safe and efficient handling.

AXS_tome⁺ Sphincterotomes (3 lumen)



Product Code	Order No.	Characteristics	Cutting Length	Tip ø Length	Accepted Guide Wire	ø Working Channel	Unit
PAP1-D7-20-35	502571	Micro-Tip ⁺ , CWT	20 mm	4,2 Fr 5 mm	.035"	2,8 mm	5
PAP1-D7-25-35	502572	Micro-Tip ⁺ , CWT	25 mm	4,2 Fr 5 mm	.035"	2,8 mm	5
PAP1-D7-30-35	502573	Micro-Tip ⁺ , CWT	30 mm	4,2 Fr 5 mm	.035"	2,8 mm	5
PAP1-D8-20-35	502574	Micro-Tip ⁺ , monofil, rotatable	20 mm	4,2 Fr 5 mm	.035"	2,8 mm	5
PAP1-D8-25-35	502575	Micro-Tip ⁺ , monofil, rotatable	25 mm	4,2 Fr 5 mm	.035"	2,8 mm	5
PAP1-D8-30-35	502576	Micro-Tip ⁺ , monofil, rotatable	30 mm	4,2 Fr 5 mm	.035"	2,8 mm	5
PAP1-P8-25-35	502577	Micro-Tip ⁺ , monofil, insulated cutting wire (ProCut), rotatable	12.5 12.5 mm proximal insulated	4,2 Fr 5 mm	.035"	2,8 mm	5
PAP1-P8-30-35	502578	Micro-Tip ⁺ , monofil, insulated cutting wire (ProCut), rotatable	15 15 mm proximal insulated	4,2 Fr 5 mm	.035"	2,8 mm	5

AXS_tome⁺ Sphincterotomes for **easyCure**® (3 lumen)



Product Code	Order No.	Characteristics	Cutting Length	Tip ø Length	Accepted Guide Wire	ø Working Channel	Unit
PAP1-Q7-20-35	502579	Micro-Tip ⁺ , CWT	20 mm	4,2 Fr 5 mm	.035"	3,8 mm	5
PAP1-Q7-25-35	502580	Micro-Tip ⁺ , CWT	25 mm	4,2 Fr 5 mm	.035"	3,8 mm	5
PAP1-Q7-30-35	502581	Micro-Tip ⁺ , CWT	30 mm	4,2 Fr 5 mm	.035"	3,8 mm	5
PAP1-S8-20-35	502582	Micro-Tip ⁺ , monofil, rotatable	20 mm	4,2 Fr 5 mm	.035"	3,8 mm	5
PAP1-S8-25-35	502583	Micro-Tip ⁺ , monofil, rotatable	25 mm	4,2 Fr 5 mm	.035"	3,8 mm	5
PAP1-S8-30-35	502584	Micro-Tip ⁺ , monofil, rotatable	30 mm	4,2 Fr 5 mm	.035"	3,8 mm	5
PAP1-T8-25-35	502585	Micro-Tip ⁺ , monofil, insulated cutting wire (ProCut), rotatable	12,5 12,5 mm proximal insulated	4,2 Fr 5 mm	.035"	3,8 mm	5
PAP1-T8-30-35	502586	Micro-Tip ⁺ , monofil, insulated cutting wire (ProCut), rotatable	15 15 mm proximal insulated	4,2 Fr 5 mm	.035"	3,8 mm	5



Sphincterotomes

Gentle probing
due to micro-tip design

Anatomical 3-D orientation

Precision control

AXS_tome

Conservative and safe approach

Our multi-stage catheter ensures stability during probing. Precision control of the catheter and reliable 11 o'clock orientation facilitate a conservative biliary approach.

Owing to a smooth guide wire transition to the micro-tip it is possible to minimise tissue irritation and reduce the risk of complications.

Preloaded **AXS_tome** Sphincterotomes (3 Lumen) with **GPS Wire**






	Product Code	Order No.	Characteristics	Cutting Length	Tip ø Length	Accepted Guide Wire	Length	ø Working Channel	Unit
	PAP1-X7-25-35-450-OL	502310	Micro-Tip, CWT, preloaded GPS Wire	25 mm	4,2 Fr 5 mm	.035" 450 cm	200 cm	2,8 mm	5
	PAP1-X8-30-35-260-OL	502286	Micro-Tip, monofil, preloaded GPS Wire , for short and long wire technique	30 mm	4,2 Fr 5 mm	.035" 260 cm	200 cm	3,8 mm	5
	PAP1-X8-30-35-450-OL	502287	Micro-Tip, monofil, preloaded GPS Wire	30 mm	4,2 Fr 5 mm	.035" 450 cm	200 cm	2,8 mm	5

Compacted Wire Technology (CWT)


Compacted Wire Technology combines the benefits of a mono-filament cutting wire with those of a braided cutting wire:

- ✓ High energy density for a clean cut
- ✓ Low thermal input due to the CWT structure
- ✓ Precision-control exploration

AXS_tome Sphincterotomes (2 Lumen)

	Product Code	Order No.	Characteristics	Cutting Length	Accepted Guide Wire	Length	ø Working Channel	Unit
	PAP1-C2-20-35-OL	501034	Micro-Tip, CWT	20 mm	.035"	200 cm	2,8 mm	5
	PAP1-C2-25-35-OL	501036	Micro-Tip, CWT	25 mm	.035"	200 cm	2,8 mm	5
	PAP1-C2-30-35-OL	501039	Micro-Tip, CWT	30 mm	.035"	200 cm	2,8 mm	5
	PAP1-C5-20-35-OL	501046	Micro-Tip, monofil	20 mm	.035"	200 cm	2,8 mm	5
	PAP1-C5-30-35-OL	501047	Micro-Tip, monofil	30 mm	.035"	200 cm	2,8 mm	5
	PAP1-D2-30-21-OL	501853	Micro-Tip, CWT, balloon enteroscopy	30 mm	.021"	260 cm	2,0 mm	2
	PAP1-A2-20-35-OL	501028	Precut, CWT	20 mm	.035"	200 cm	2,8 mm	5
	PAP1-K2-06-35-OL	501060	Needle	6 mm	.035"	200 cm	2,8 mm	5
	PAP1-M2-20-35-OL	501840	B II, CWT	20 mm	.035"	200 cm	2,8 mm	2
	PAP1-N2-20-21-OL	501902	B II, CWT, balloon enteroscopy	20 mm	.021"	260 cm	2,0 mm	2

Anatomical 3-D orientation
Based on systematic 11 o'clock orientation
the AXS_tome sphincterotomes also follow
the left frontal anatomy of the
distal choledochus.



AXS_tome Sphincterotomes for **easycure**® (Short and long wire | 3 Lumen)




	Product Code	Order No.	Characteristics	Cutting Length	Tip ø Length	Accepted Guide Wire	Length	ø Working Channel	Unit
	PAP1-K8-06-E35-OL	501920	Needle	6 mm	-	.035" ≥ 200 cm	200 cm	3,8 mm	5
	PAP1-X8-30-35-260-OL	502286	Micro-Tip, monofil, preloaded GPS Wire	30 mm	4,2 Fr 5 mm	.035" 260 cm	200 cm	3,8 mm	5

easycure®: The platform
increasing your options

Where other systems tend to restrict you, easycure offers extra options for ERCP. No matter whether it is a question of short or long-wire systems, the decision is yours. All the instruments are designed to allow direct guide wire securement, either on the biopsy inlet channel or using an Albarran lever.





	Product Code	Order No.	Characteristics	Unit
	SPE1-Y7	501925	Guide wire fixation (instrument fork), for Olympus and Fujifilm (from G5 series on)	5

ERCP cannulas



	Product Code	Order No.	Characteristics	Length	Accepted Guide Wire	ø Working Channel	Unit
	CAN1-A3-18-200-35	500337	Standard, 1 lumen, radiopaque marking	200 cm	.035"	2,0 mm	10
	CAN1-A5-18-200-35	500339	Filiiform, 1 lumen, radiopaque marking	200 cm	.035"	2,0 mm	10
	CAN1-A6-18-200-35	502332	Filiiform, guide wire port, 1 lumen, radiopaque marking, tapered (2,3 1,8 mm)	200 cm	.035"	2,8 mm	10
	CAN1-B3-18-200-35	500343	Filiiform, 1 lumen, radiopaque marking	200 cm	.035"	2,0 mm	10
	CAN1-B3-18-260-35	501856	Filiiform, 1 lumen, radiopaque marking, balloon enteroscopy	260 cm	.035"	2,0 mm	2
	CAN1-B4-18-200-35	502008	Filiiform, 2 lumen, radiopaque marking, tapered (2,3 1,8 mm)	200 cm	.035"	2,8 mm	10
	CAN1-B4-18-200-E35	502007	Filiiform, 2 lumen, radiopaque marking, for short and long wire technique, tapered (2,3 1,8 mm)	200 cm	.035" ≥ 200 cm	3,2 mm	2
	CAN1-D5-18-200-35	500348	Bottle neck, 1 lumen, metal tip	200 cm	.035"	2,0 mm	10



Guide wires

Precise exploration

Stable instrument guidance

Smooth transition between tip and coating

GPS Wire

Arrive safely

The atraumatic tip of the GI Probing Scout ensures conservative exploration. The sturdy sheath reliably guides therapy instruments to their destination. A smooth transition

between the hydrophilic tip and the coating enhances the precision control of your instruments. In addition, handling is facilitated by the newly designed insertion aid.

GPS Wire Guide wires with hydrophilic tip



	Product Code	Order No.	Characteristics	Diameter	Length	Unit
	WIR1-S6-35-260	502517	Straight, Nitinol, insulated, green/white	.035"	260 cm	5
	WIR1-S6-35-450	502515	Straight, Nitinol, insulated, green/white	.035"	450 cm	5
	WIR1-J6-35-260	502518	J-tip, Nitinol, insulated, green/white	.035"	260 cm	5
	WIR1-J6-35-450	502516	J-tip, Nitinol, insulated, green/white	.035"	450 cm	5


Guide wires



	Product Code	Order No.	Characteristics	Diameter	Length	Unit
	WIR1-E1-72-260	501780	Straight, steel, flexible tip, Eder-Püstow	.072"	260 cm	2
	WIR1-F1-18-400	501781	Straight, Nitinol, coated	.018"	400 cm	2
	WIR1-F1-21-400	501782	Straight, Nitinol, coated	.021"	400 cm	2
	WIR1-F1-21-500	501911	Straight, Nitinol, coated	.021"	500 cm	2
	WIR1-F5-25-400	501789	Straight, Nitinol, insulated, black/white	.025"	400 cm	2
	WIR1-F5-35-400	501792	Straight, Nitinol, insulated, black/white	.035"	400 cm	2
	WIR1-F7-35-260	501837	Straight, Nitinol, coated, extra stiff	.035"	260 cm	2
	WIR1-F7-35-400	501838	Straight, Nitinol, coated, extra stiff	.035"	400 cm	2

1-9

Stone extraction balloons



Product Code	Order No.	Characteristics	ø Balloon	Length	Accepted Guide Wire	ø Working Channel	Unit
BAL1-B2-16-75-200-35	500005	Distal radiopaque marking, 2 lumen, tapered (7 5 Fr)	16 mm	200 cm	.035"	2,8 mm	5
BAL1-C2-16-07-200-35	500010	Distal radiopaque marking, 2 lumen	16 mm	200 cm	.035"	2,8 mm	5
BAL1-C3-18-07-200-35	500011	Distal radiopaque marking, 3 lumen	16 mm	200 cm	.035"	2,8 mm	5
BAL1-C3-18-07-200E35	501916	Distal radiopaque marking, tapered, 3 lumen, for short and long wire technique	16 mm	200 cm	.035" ≥ 200 cm	3,2 mm	2



Maximum expansion force

Precision control by rotary motion

Infinite adjustment of basket size

twist'n'CATCH

Collects and retains reliably


The impressive expansion force and precise control of the basket allow swift capture of stones. Infinite adjustment of basket

diameter permits highly flexible use. As a result, gallstones can be retrieved reliably.


twist'n'CATCH Stone extraction baskets



Product Code	Order No.	Characteristics	ø Basket	Length	Accepted Guide Wire	ø Working Channel	Unit
BAS1-K2-25-23-200	500027	Spherical, 130° rotation	35 mm	200 cm	-	2,8 mm	5
BAS1-K3-35-28-200	502285	Spherical, 130° rotation, for short and long wire technique	35 mm	200 cm	.035" ≥ 200 cm	3,2 mm	5



CAUTION: FUJIFILM medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of FUJIFILM medwork for such cases is excluded.



For guide wires $\geq 200^*$



Reliable expansion even in challenging applications

*Figure: BAS1-K3-35-28-200




The right tool for every situation

No matter how large the ,rock‘ may be, Rockstar is fearless: with the new stone remover for lithotripsy you are ideally equipped for retrieving gallstones.

The guide wire lumen enables you to reach the location of the stone reliably. The innovative operating concept allows intuitive pulverisation of concretions if necessary.

ROCKSTAR Lithotriptor stone extraction basket




	Product Code	Order No.	Characteristics	ø Basket	Length	Accepted Guide Wire	ø Working Channel	Unit
	BAS3-K2-35-35-200	502486	twist'n'CATCH-basket, integrated lithotripsy function	35 mm	200 cm	.035"	4,2 mm	2

CAUTION: medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Stone extraction baskets



	Product Code	Order No.	Characteristics	ø Basket	Length	ø Working Channel	Unit
	BAS1-A2-20-23-200	500018	Hexagonal	20 mm	200 cm	2,8 mm	5
	BAS1-A2-25-23-200	500020	Hexagonal	25 mm	200 cm	2,8 mm	5
	BAS1-A2-30-23-200	500022	Hexagonal	30 mm	200 cm	2,8 mm	5
	BAS1-A2-35-23-200	500024	Hexagonal	35 mm	200 cm	2,8 mm	5
	BAS1-C1-20-18-260	500025	Half twisted, balloon enteroscopy	20 mm	260 cm	2,0 mm	2
	BAS1-C2-25-23-260	500026	Half twisted, balloon enteroscopy	25 mm	260 cm	2,8 mm	2

CAUTION: medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

ERCP

Lithotripter baskets



	Product Code	Order No.	Characteristics	ø Basket	Length	ø Working Channel	Unit
	twist'n'CATCH	500826	Spherical, 130° rotation	35 mm	200 400 cm	2,8 mm	2
	LIT1-K2-25-23-200						
	LIT1-C3-20-23-200	500822	Half twisted	20 mm	200 400 cm	2,8 mm	2
	LIT1-C3-25-23-200	500823	Half twisted	25 mm	200 400 cm	2,8 mm	2
	LIT1-C3-30-23-200	500824	Half twisted	30 mm	200 400 cm	2,8 mm	2
	LIT1-C3-35-23-200	500825	Half twisted	35 mm	200 400 cm	2,8 mm	2

CAUTION: medwork products of the LIT series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Lithotripter springs



	Product Code	Order No.	Characteristics	Length	ø Working Channel	Unit
	LIT1-G4-26-185	502304	Spiral spring, transendoscopic, bile duct	185 cm	2,8 mm	2
	LIT1-N4-38-090	502305	Spiral spring, emergency, bile duct	90 cm	-	2

CAUTION: medwork products of the LIT series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Handle (Reusable)

134°C

	Product Code	Order No.	Characteristics	Unit
	LIT2-D-HAN	500881	Handle, pull-through technique	1

CAUTION: medwork products of the LIT series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.



#flowing

Infinitely adjustable guide catheter length

Precise placement of prosthesis
through accurate force transmission

The placement system: One-for-All

KENA

Harmonically combined

Inspired by the music, tuned to the finest detail: the new Kena series allows simple and precise placement of the prosthesis as a result of the excellent interplay between bending strength and flexibility.

Even multi-stenting with different prosthetic diameters (8.5 and 10 Fr), including repositioning, is extremely efficient with the One-for-All placement system.

KENA Biliary stents (2 flaps, centre bend)



	Product Code	Order No.	Color	Diameter	Length Between Flaps	Total Length	Unit
	PR01-C2-07-05	502362	Grey	7 Fr	5 cm	6,5 cm	1
	PR01-C2-07-07	502368	Grey	7 Fr	7 cm	8,5 cm	1
	PR01-C2-07-09	502366	Grey	7 Fr	9 cm	10,5 cm	1
	PR01-C2-07-11	502367	Grey	7 Fr	11 cm	12,5 cm	1
	PR01-C2-07-13	502369	Grey	7 Fr	13 cm	14,5 cm	1
	PR01-C2-07-15	502370	Grey	7 Fr	15 cm	16,5 cm	1
	PR01-C2-07-17	502371	Grey	7 Fr	17 cm	18,5 cm	1
	PR01-C2-07-19	502372	Grey	7 Fr	19 cm	20,5 cm	1
	PR01-C2-08-05	502373	Green	8.5 Fr	5 cm	6,5 cm	1
	PR01-C2-08-07	502374	Green	8.5 Fr	7 cm	8,5 cm	1
	PR01-C2-08-09	502375	Green	8.5 Fr	9 cm	10,5 cm	1
	PR01-C2-08-11	502376	Green	8.5 Fr	11 cm	12,5 cm	1
	PR01-C2-08-13	502377	Green	8.5 Fr	13 cm	14,5 cm	1
	PR01-C2-08-15	502378	Green	8.5 Fr	15 cm	16,5 cm	1
	PR01-C2-08-17	502379	Green	8.5 Fr	17 cm	18,5 cm	1
	PR01-C2-08-19	502380	Green	8.5 Fr	19 cm	20,5 cm	1
	PR01-C2-10-05	502381	Blue	10 Fr	5 cm	6,5 cm	1
	PR01-C2-10-06	502419	Blue	10 Fr	6 cm	7,5 cm	1
	PR01-C2-10-07	502382	Blue	10 Fr	7 cm	8,5 cm	1
	PR01-C2-10-08	502420	Blue	10 Fr	8 cm	9,5 cm	1
	PR01-C2-10-09	502383	Blue	10 Fr	9 cm	10,5 cm	1
	PR01-C2-10-10	502421	Blue	10 Fr	10 cm	11,5 cm	1
	PR01-C2-10-11	502384	Blue	10 Fr	11 cm	12,5 cm	1
	PR01-C2-10-12	502422	Blue	10 Fr	12 cm	13,5 cm	1
	PR01-C2-10-13	502385	Blue	10 Fr	13 cm	14,5 cm	1
	PR01-C2-10-14	502423	Blue	10 Fr	14 cm	15,5 cm	1
	PR01-C2-10-15	502386	Blue	10 Fr	15 cm	16,5 cm	1
	PR01-C2-10-17	502387	Blue	10 Fr	17 cm	18,5 cm	1
	PR01-C2-10-19	502388	Blue	10 Fr	19 cm	20,5 cm	1

KENA Biliary stents (2 flaps, proximal bend)






Product Code	Order No.	Color	Diameter	Length Between Flaps	Total Length	Unit
PR01-T2-07-05	502392	Grey	7 Fr	5 cm	6,5 cm	1
PR01-T2-07-07	502395	Grey	7 Fr	7 cm	8,5 cm	1
PR01-T2-07-09	502393	Grey	7 Fr	9 cm	10,5 cm	1
PR01-T2-07-11	502394	Grey	7 Fr	11 cm	12,5 cm	1
PR01-T2-07-13	502396	Grey	7 Fr	13 cm	14,5 cm	1
PR01-T2-07-15	502397	Grey	7 Fr	15 cm	16,5 cm	1
PR01-T2-07-17	502398	Grey	7 Fr	17 cm	18,5 cm	1
PR01-T2-07-19	502399	Grey	7 Fr	19 cm	20,5 cm	1
PR01-T2-08-05	502400	Green	8.5 Fr	5 cm	6,5 cm	1
PR01-T2-08-07	502401	Green	8.5 Fr	7 cm	8,5 cm	1
PR01-T2-08-09	502402	Green	8.5 Fr	9 cm	10,5 cm	1
PR01-T2-08-11	502403	Green	8.5 Fr	11 cm	12,5 cm	1
PR01-T2-08-13	502404	Green	8.5 Fr	13 cm	14,5 cm	1
PR01-T2-08-15	502405	Green	8.5 Fr	15 cm	16,5 cm	1
PR01-T2-08-17	502406	Green	8.5 Fr	17 cm	18,5 cm	1
PR01-T2-08-19	502407	Green	8.5 Fr	19 cm	20,5 cm	1
PR01-T2-10-05	502408	Blue	10 Fr	5 cm	6,5 cm	1
PR01-T2-10-06	502424	Blue	10 Fr	6 cm	7,5 cm	1
PR01-T2-10-07	502409	Blue	10 Fr	7 cm	8,5 cm	1
PR01-T2-10-08	502425	Blue	10 Fr	8 cm	9,5 cm	1
PR01-T2-10-09	502410	Blue	10 Fr	9 cm	10,5 cm	1
PR01-T2-10-10	502426	Blue	10 Fr	10 cm	11,5 cm	1
PR01-T2-10-11	502411	Blue	10 Fr	11 cm	12,5 cm	1
PR01-T2-10-12	502427	Blue	10 Fr	12 cm	13,5 cm	1
PR01-T2-10-13	502412	Blue	10 Fr	13 cm	14,5 cm	1
PR01-T2-10-14	502428	Blue	10 Fr	14 cm	15,5 cm	1
PR01-T2-10-15	502413	Blue	10 Fr	15 cm	16,5 cm	1
PR01-T2-10-17	502414	Blue	10 Fr	17 cm	18,5 cm	1
PR01-T2-10-19	502415	Blue	10 Fr	19 cm	20,5 cm	1

Preloaded **KENA** Biliary stents (2 flaps, centre bend)

For long wire technique




Product Code	Order No.	Color	Diameter	Length Between Flaps	Total Length	ø Working Channel	Unit
PR03-C2-08-05	502435	Green	8.5 Fr	5 cm	6,5 cm	4,2 mm	1
PR03-C2-08-07	502438	Green	8.5 Fr	7 cm	8,5 cm	4,2 mm	1
PR03-C2-08-09	502439	Green	8.5 Fr	9 cm	10,5 cm	4,2 mm	1
PR03-C2-08-11	502440	Green	8.5 Fr	11 cm	12,5 cm	4,2 mm	1
PR03-C2-08-13	502441	Green	8.5 Fr	13 cm	14,5 cm	4,2 mm	1
PR03-C2-08-15	502442	Green	8.5 Fr	15 cm	16,5 cm	4,2 mm	1
PR03-C2-08-17	502443	Green	8.5 Fr	17 cm	18,5 cm	4,2 mm	1
PR03-C2-08-19	502502	Green	8.5 Fr	19 cm	20,5 cm	4,2 mm	1
PR03-C2-10-05	502444	Blue	10 Fr	5 cm	6,5 cm	4,2 mm	1
PR03-C2-10-06	502445	Blue	10 Fr	6 cm	7,5 cm	4,2 mm	1
PR03-C2-10-07	502446	Blue	10 Fr	7 cm	8,5 cm	4,2 mm	1
PR03-C2-10-08	502447	Blue	10 Fr	8 cm	9,5 cm	4,2 mm	1
PR03-C2-10-09	502448	Blue	10 Fr	9 cm	10,5 cm	4,2 mm	1
PR03-C2-10-10	502449	Blue	10 Fr	10 cm	11,5 cm	4,2 mm	1
PR03-C2-10-11	502450	Blue	10 Fr	11 cm	12,5 cm	4,2 mm	1
PR03-C2-10-12	502451	Blue	10 Fr	12 cm	13,5 cm	4,2 mm	1
PR03-C2-10-13	502452	Blue	10 Fr	13 cm	14,5 cm	4,2 mm	1
PR03-C2-10-14	502453	Blue	10 Fr	14 cm	15,5 cm	4,2 mm	1
PR03-C2-10-15	502454	Blue	10 Fr	15 cm	16,5 cm	4,2 mm	1
PR03-C2-10-17	502455	Blue	10 Fr	17 cm	18,5 cm	4,2 mm	1
PR03-C2-10-19	502503	Blue	10 Fr	19 cm	20,5 cm	4,2 mm	1
PR03-C2-10-21	502504	Blue	10 Fr	21 cm	22,5 cm	4,2 mm	1
PR03-C2-10-23	502505	Blue	10 Fr	23 cm	24,5 cm	4,2 mm	1
PR03-C2-10-25	502506	Blue	10 Fr	25 cm	26,5 cm	4,2 mm	1

Preloaded **KENA** Biliary stents (2 flaps, proximal bend)

For long wire technique




Product Code	Order No.	Color	Diameter	Length Between Flaps	Total Length	ø Working Channel	Unit
PR03-T2-08-05	502436	Green	8.5 Fr	5 cm	6,5 cm	4,2 mm	1
PR03-T2-08-07	502456	Green	8.5 Fr	7 cm	8,5 cm	4,2 mm	1
PR03-T2-08-09	502457	Green	8.5 Fr	9 cm	10,5 cm	4,2 mm	1
PR03-T2-08-11	502458	Green	8.5 Fr	11 cm	12,5 cm	4,2 mm	1
PR03-T2-08-13	502459	Green	8.5 Fr	13 cm	14,5 cm	4,2 mm	1
PR03-T2-08-15	502460	Green	8.5 Fr	15 cm	16,5 cm	4,2 mm	1
PR03-T2-08-17	502461	Green	8.5 Fr	17 cm	18,5 cm	4,2 mm	1
PR03-T2-08-19	502507	Green	8.5 Fr	19 cm	20,5 cm	4,2 mm	1
PR03-T2-10-05	502462	Blue	10 Fr	5 cm	6,5 cm	4,2 mm	1
PR03-T2-10-06	502463	Blue	10 Fr	6 cm	7,5 cm	4,2 mm	1
PR03-T2-10-07	502464	Blue	10 Fr	7 cm	8,5 cm	4,2 mm	1
PR03-T2-10-08	502465	Blue	10 Fr	8 cm	9,5 cm	4,2 mm	1
PR03-T2-10-09	502466	Blue	10 Fr	9 cm	10,5 cm	4,2 mm	1
PR03-T2-10-10	502467	Blue	10 Fr	10 cm	11,5 cm	4,2 mm	1
PR03-T2-10-11	502468	Blue	10 Fr	11 cm	12,5 cm	4,2 mm	1
PR03-T2-10-12	502469	Blue	10 Fr	12 cm	13,5 cm	4,2 mm	1
PR03-T2-10-13	502470	Blue	10 Fr	13 cm	14,5 cm	4,2 mm	1
PR03-T2-10-14	502471	Blue	10 Fr	14 cm	15,5 cm	4,2 mm	1
PR03-T2-10-15	502472	Blue	10 Fr	15 cm	16,5 cm	4,2 mm	1
PR03-T2-10-17	502473	Blue	10 Fr	17 cm	18,5 cm	4,2 mm	1
PR03-T2-10-19	502508	Blue	10 Fr	19 cm	20,5 cm	4,2 mm	1
PR03-T2-10-21	502509	Blue	10 Fr	21 cm	22,5 cm	4,2 mm	1
PR03-T2-10-23	502510	Blue	10 Fr	23 cm	24,5 cm	4,2 mm	1
PR03-T2-10-25	502511	Blue	10 Fr	25 cm	26,5 cm	4,2 mm	1




Tips tapered on both sides

The unique bilateral tip facilitates placement and minimizes the perforation risk. Alternatively, the prosthesis can be laid across a bifurcation with the bend first and thus be anchored in place. This ensures a better hold and reduces the risk of migration.



Preloaded **KENA** Biliary stents (2 flaps, centre bend) **easyCure**

For short wire technique




Product Code	Order No.	Color	Diameter	Length Between Flaps	Total Length	ø Working Channel	Unit
PR03-C3-08-05	502652	Green	8.5 Fr	5 cm	6,5 cm	4,2 mm	1
PR03-C3-08-07	502653	Green	8.5 Fr	7 cm	8,5 cm	4,2 mm	1
PR03-C3-08-09	502654	Green	8.5 Fr	9 cm	10,5 cm	4,2 mm	1
PR03-C3-08-11	502655	Green	8.5 Fr	11 cm	12,5 cm	4,2 mm	1
PR03-C3-08-13	502656	Green	8.5 Fr	13 cm	14,5 cm	4,2 mm	1
PR03-C3-08-15	502657	Green	8.5 Fr	15 cm	16,5 cm	4,2 mm	1
PR03-C3-08-17	502658	Green	8.5 Fr	17 cm	18,5 cm	4,2 mm	1
PR03-C3-08-19	502659	Green	8.5 Fr	19 cm	20,5 cm	4,2 mm	1
PR03-C3-10-05	502660	Blue	10 Fr	5 cm	6,5 cm	4,2 mm	1
PR03-C3-10-06	502661	Blue	10 Fr	6 cm	7,5 cm	4,2 mm	1
PR03-C3-10-07	502662	Blue	10 Fr	7 cm	8,5 cm	4,2 mm	1
PR03-C3-10-08	502663	Blue	10 Fr	8 cm	9,5 cm	4,2 mm	1
PR03-C3-10-09	502664	Blue	10 Fr	9 cm	10,5 cm	4,2 mm	1
PR03-C3-10-10	502665	Blue	10 Fr	10 cm	11,5 cm	4,2 mm	1
PR03-C3-10-11	502666	Blue	10 Fr	11 cm	12,5 cm	4,2 mm	1
PR03-C3-10-12	502667	Blue	10 Fr	12 cm	13,5 cm	4,2 mm	1
PR03-C3-10-13	502668	Blue	10 Fr	13 cm	14,5 cm	4,2 mm	1
PR03-C3-10-14	502669	Blue	10 Fr	14 cm	15,5 cm	4,2 mm	1
PR03-C3-10-15	502670	Blue	10 Fr	15 cm	16,5 cm	4,2 mm	1
PR03-C3-10-17	502671	Blue	10 Fr	17 cm	18,5 cm	4,2 mm	1
PR03-C3-10-19	502672	Blue	10 Fr	19 cm	20,5 cm	4,2 mm	1
PR03-C3-10-21	502673	Blue	10 Fr	21 cm	22,5 cm	4,2 mm	1
PR03-C3-10-23	502674	Blue	10 Fr	23 cm	24,5 cm	4,2 mm	1
PR03-C3-10-25	502675	Blue	10 Fr	25 cm	26,5 cm	4,2 mm	1

Preloaded **KENA** Biliary stents (2 flaps, proximal bend) *easyCure*[®]

For short wire technique




	Product Code	Order No.	Color	Diameter	Length Between Flaps	Total Length	ø Working Channel	Unit
	PR03-T3-08-05	502676	Green	8.5 Fr	5 cm	6,5 cm	4,2 mm	1
	PR03-T3-08-07	502677	Green	8.5 Fr	7 cm	8,5 cm	4,2 mm	1
	PR03-T3-08-09	502678	Green	8.5 Fr	9 cm	10,5 cm	4,2 mm	1
	PR03-T3-08-11	502679	Green	8.5 Fr	11 cm	12,5 cm	4,2 mm	1
	PR03-T3-08-13	502680	Green	8.5 Fr	13 cm	14,5 cm	4,2 mm	1
	PR03-T3-08-15	502681	Green	8.5 Fr	15 cm	16,5 cm	4,2 mm	1
	PR03-T3-08-17	502682	Green	8.5 Fr	17 cm	18,5 cm	4,2 mm	1
	PR03-T3-08-19	502683	Green	8.5 Fr	19 cm	20,5 cm	4,2 mm	1
	PR03-T3-10-05	502684	Blue	10 Fr	5 cm	6,5 cm	4,2 mm	1
	PR03-T3-10-06	502685	Blue	10 Fr	6 cm	7,5 cm	4,2 mm	1
	PR03-T3-10-07	502686	Blue	10 Fr	7 cm	8,5 cm	4,2 mm	1
	PR03-T3-10-08	502687	Blue	10 Fr	8 cm	9,5 cm	4,2 mm	1
	PR03-T3-10-09	502688	Blue	10 Fr	9 cm	10,5 cm	4,2 mm	1
	PR03-T3-10-10	502689	Blue	10 Fr	10 cm	11,5 cm	4,2 mm	1
	PR03-T3-10-11	502690	Blue	10 Fr	11 cm	12,5 cm	4,2 mm	1
	PR03-T3-10-12	502691	Blue	10 Fr	12 cm	13,5 cm	4,2 mm	1
	PR03-T3-10-13	502692	Blue	10 Fr	13 cm	14,5 cm	4,2 mm	1
	PR03-T3-10-14	502693	Blue	10 Fr	14 cm	15,5 cm	4,2 mm	1
	PR03-T3-10-15	502694	Blue	10 Fr	15 cm	16,5 cm	4,2 mm	1
	PR03-T3-10-17	502695	Blue	10 Fr	17 cm	18,5 cm	4,2 mm	1
	PR03-T3-10-19	502696	Blue	10 Fr	19 cm	20,5 cm	4,2 mm	1
	PR03-T3-10-21	502697	Blue	10 Fr	21 cm	22,5 cm	4,2 mm	1
	PR03-T3-10-23	502698	Blue	10 Fr	23 cm	24,5 cm	4,2 mm	1
	PR03-T3-10-25	502699	Blue	10 Fr	25 cm	26,5 cm	4,2 mm	1

KENA Pusher

For short and long wire technique




	Product Code	Order No.	Characteristics	Diameter	Length	Accepted Guide Wire	ø Working Channel	Unit
	PR01-P2-07-180	502494	Pusher, white, for 7 Fr biliary stents	7 Fr	180 cm	.035"	2,8 mm	2
	PR01-P1-09-180	502492	Pusher, white, for 8.5 and 10 Fr biliary stents	9 Fr	180 cm	.035"	4,2 mm	2

KENA Guide catheter and insertion kit

For long wire technique




	Product Code	Order No.	Characteristics	Diameter	Length	Accepted Guide Wire	ø Working Channel	Unit
	PR01-K1-05-220-35	502491	Guide catheter, grey, for 8.5 and 10 Fr biliary stents	5 Fr	220 cm	.035"	4,2 mm	2
	PR01-S1-05-220-35	502493	Insertion kit for 8.5 and 10 Fr biliary stents, LL-connection, distal radiopaque marking	9 Fr	220 cm	.035"	4,2 mm	2

KENA Guide catheter and insertion kit

For short wire technique



	Product Code	Order No.	Characteristics	Diameter	Length	Accepted Guide Wire	ø Working Channel	Unit
	PR01-K3-05-220-35	502650	Guide catheter, grey, for 8.5 and 10 Fr biliary stents	5 Fr	220 cm	.035"	4,2 mm	2
	PR01-S1-09-220-35	502651	Insertion kit for 8.5 and 10 Fr biliary stents, LL-connection, distal radiopaque marking	9 Fr	220 cm	.035"	4,2 mm	2


Dilation balloons

Working Length 230 cm



Product Code	Order No.	Characteristics	ø Balloon	Length Balloon	Accepted Guide Wire	ø Working Channel	Unit
DIL1-A1-06-40-23-35	502268	2 radiopaque markings, tapered catheter (2,3 1,8 mm), for short and long wire technique	6 mm	40 mm	.035" ≥ 200 cm	2,8 mm	2
DIL1-A1-08-40-23-35	502269	2 radiopaque markings, tapered catheter (2,3 1,8 mm), for short and long wire technique	8 mm	40 mm	.035" ≥ 200 cm	2,8 mm	2
DIL1-A1-10-40-23-35	502270	2 radiopaque markings, tapered catheter (2,3 1,8 mm), for short and long wire technique	10 mm	40 mm	.035" ≥ 200 cm	2,8 mm	2
DIL1-A1-12-60-23-35	502271	2 radiopaque markings, for short and long wire technique	12 mm	60 mm	.035" ≥ 200 cm	2,8 mm	2
DIL1-A1-14-60-23-35	502272	2 radiopaque markings, for short and long wire technique	14 mm	60 mm	.035" ≥ 200 cm	2,8 mm	2
DIL1-A1-16-60-23-35	502273	2 radiopaque markings, for short and long wire technique	16 mm	60 mm	.035" ≥ 200 cm	2,8 mm	2
DIL1-A1-18-60-23-35	502274	2 radiopaque markings, for short and long wire technique	18 mm	60 mm	.035" ≥ 200 cm	2,8 mm	2
DIL1-A1-20-60-23-35	502275	2 radiopaque markings, for short and long wire technique	20 mm	60 mm	.035" ≥ 200 cm	3,2 mm	2

Insufflation device



Product Code	Order No.	Characteristics	Unit
DIL1-B1-60	502278	60 ml volume, 12 atm maximum pressure	1

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231	Калининград (4012)72-03-81	Омск (3812)21-46-40	Сыктывкар (8212)25-95-17
Ангарск (3955)60-70-56	Калуга (4842)92-23-67	Орел (4862)44-53-42	Тамбов (4752)50-40-97
Архангельск (8182)63-90-72	Кемерово (3842)65-04-62	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Астрахань (8512)99-46-04	Киров (8332)68-02-04	Пенза (8412)22-31-16	Тольятти (8482)63-91-07
Барнаул (3852)73-04-60	Коломна (4966)23-41-49	Петрозаводск (8142)55-98-37	Томск (3822)98-41-53
Белгород (4722)40-23-64	Кострома (4942)77-07-48	Псков (8112)59-10-37	Тула (4872)33-79-87
Благовещенск (4162)22-76-07	Краснодар (861)203-40-90	Пермь (342)205-81-47	Тюмень (3452)66-21-18
Брянск (4832)59-03-52	Красноярск (391)204-63-61	Ростов-на-Дону (863)308-18-15	Ульяновск (8422)24-23-59
Владивосток (423)249-28-31	Курск (4712)77-13-04	Рязань (4912)46-61-64	Улан-Удэ (3012)59-97-51
Владикавказ (8672)28-90-48	Курган (3522)50-90-47	Самара (846)206-03-16	Уфа (347)229-48-12
Владимир (4922)49-43-18	Липецк (4742)52-20-81	Саранск (8342)22-96-24	Хабаровск (4212)92-98-04
Волгоград (844)278-03-48	Магнитогорск (3519)55-03-13	Санкт-Петербург (812)309-46-40	Чебоксары (8352)28-53-07
Вологда (8172)26-41-59	Москва (495)268-04-70	Саратов (845)249-38-78	Челябинск (351)202-03-61
Воронеж (473)204-51-73	Мурманск (8152)59-64-93	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Екатеринбург (343)384-55-89	Набережные Челны (8552)20-53-41	Симферополь (3652)67-13-56	Чита (3022)38-34-83
Иваново (4932)77-34-06	Нижний Новгород (831)429-08-12	Смоленск (4812)29-41-54	Якутск (4112)23-90-97
Ижевск (3412)26-03-58	Новокузнецк (3843)20-46-81	Сочи (862)225-72-31	Ярославль (4852)69-52-93
Иркутск (395)279-98-46	Ноябрьск (3496)41-32-12	Ставрополь (8652)20-65-13	
Казань (843)206-01-48	Новосибирск (383)227-86-73	Сургут (3462)77-98-35	
Россия +7(495)268-04-70	Киргизия +996(312)-96-26-47	Казахстан +7(7172)727-132	