

Терапевтические устройства и аксессуары

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48

Россия +7(495)268-04-70

Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81
Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73

Киргизия +996(312)-96-26-47

Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35


Казахстан +7(7172)727-132

Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Available as of January 2020



AXS_tome⁺ Sphincterotomes (3 lumen)




Product Code	Order No.	Characteristics	Cutting Length	Tip ø Length	Accepted Guide Wire	ø Working Channel	Unit
PAP1-D7-20-35	502571	Micro-Tip ⁺ , CWT	20 mm	4.2 Fr 5 mm	.035"	2.8 mm	5
PAP1-D7-25-35	502572	Micro-Tip ⁺ , CWT	25 mm	4.2 Fr 5 mm	.035"	2.8 mm	5
PAP1-D7-30-35	502573	Micro-Tip ⁺ , CWT	30 mm	4.2 Fr 5 mm	.035"	2.8 mm	5
PAP1-D8-20-35	502574	Micro-Tip ⁺ , monofil, rotatable	20 mm	4.2 Fr 5 mm	.035"	2.8 mm	5
PAP1-D8-25-35	502575	Micro-Tip ⁺ , monofil, rotatable	25 mm	4.2 Fr 5 mm	.035"	2.8 mm	5
PAP1-D8-30-35	502576	Micro-Tip ⁺ , monofil, rotatable	30 mm	4.2 Fr 5 mm	.035"	2.8 mm	5
PAP1-P8-25-35	502577	Micro-Tip ⁺ , monofil, insulated cutting wire (ProCut), rotatable	12.5 12.5 mm proximal insulated	4.2 Fr 5 mm	.035"	2.8 mm	5
PAP1-P8-30-35	502578	Micro-Tip ⁺ , monofil, insulated cutting wire (ProCut), rotatable	15 15 mm proximal insulated	4.2 Fr 5 mm	.035"	2.8 mm	5

Available as of January 2020




AXS_tome⁺ Sphincterotomes for **easyCure**® (3 lumen)



Product Code	Order No.	Characteristics	Cutting Length	Tip ø Length	Accepted Guide Wire	ø Working Channel	Unit
PAP1-Q7-20-35	502579	Micro-Tip ⁺ , CWT	20 mm	4.2 Fr 5 mm	.035"	3.8 mm	5
PAP1-Q7-25-35	502580	Micro-Tip ⁺ , CWT	25 mm	4.2 Fr 5 mm	.035"	3.8 mm	5
PAP1-Q7-30-35	502581	Micro-Tip ⁺ , CWT	30 mm	4.2 Fr 5 mm	.035"	3.8 mm	5
PAP1-S8-20-35	502582	Micro-Tip ⁺ , monofil, rotatable	20 mm	4.2 Fr 5 mm	.035"	3.8 mm	5
PAP1-S8-25-35	502583	Micro-Tip ⁺ , monofil, rotatable	25 mm	4.2 Fr 5 mm	.035"	3.8 mm	5
PAP1-S8-30-35	502584	Micro-Tip ⁺ , monofil, rotatable	30 mm	4.2 Fr 5 mm	.035"	3.8 mm	5
PAP1-T8-25-35	502585	Micro-Tip ⁺ , monofil, insulated cutting wire (ProCut), rotatable	12.5 12.5 mm proximal insulated	4.2 Fr 5 mm	.035"	3.8 mm	5
PAP1-T8-30-35	502586	Micro-Tip ⁺ , monofil, insulated cutting wire (ProCut), rotatable	15 15 mm proximal insulated	4.2 Fr 5 mm	.035"	3.8 mm	5



Preloaded **AXS_tome** Sphincterotomes (3 lumen) with gSlider



	Product Code	Order No.	Characteristics	Cutting Length	gSlider	ø Working Channel	Unit
	PAP1-X7-25-35-450-OL	502310	Micro-Tip, CWT, preloaded gSlider	25 mm	.035" 450 cm	2.8 mm	5
	PAP1-X8-30-35-260-OL	502286	Micro-Tip, monofil, preloaded gSlider, for short and long wire technique	30 mm	.035" 260 cm	3.8 mm	5
	PAP1-X8-30-35-450-OL	502287	Micro-Tip, monofil, preloaded gSlider	30 mm	.035" 450 cm	2.8 mm	5

AXS_tome Sphincterotomes (3 lumen)




	Product Code	Order No.	Characteristics	Cutting Length	Accepted Guide Wire	ø Working Channel	Unit
	PAP1-C7-20-35-OL	501048	Micro-Tip, CWT	20 mm	.035"	2.8 mm	5
	PAP1-C7-25-35-OL	502290	Micro-Tip, CWT	25 mm	.035"	2.8 mm	5
	PAP1-C7-30-35-OL	501050	Micro-Tip, CWT	30 mm	.035"	2.8 mm	5
	PAP1-C8-20-35-OL	501051	Micro-Tip, monofil	20 mm	.035"	2.8 mm	5
	PAP1-C8-25-35-OL	502291	Micro-Tip, monofil	25 mm	.035"	2.8 mm	5
	PAP1-C8-30-35-OL	501052	Micro-Tip, monofil	30 mm	.035"	2.8 mm	5
	PAP1-F8-25-35-OL	502501	Micro-Tip, monofil, insulated cutting wire (ProCut)	12.5 12.5 mm proximal insulated	.035"	2.8 mm	5
	PAP1-F8-30-35-OL	501993	Micro-Tip, monofil, insulated cutting wire (ProCut)	15 15 mm proximal insulated	.035"	2.8 mm	5

Compacted Wire Technologie (CWT)


Compacted Wire Technology combines the benefits of a monofilament cutting wire with those of a braided cutting wire:

- ✓ High energy density for a clean cut
- ✓ Low thermal input due to the CWT structure
- ✓ Precision-control probing

AXS_tome Sphincterotomes (2 lumen)




Product Code	Order No.	Characteristics	Cutting Length	Accepted Guide Wire	ø Working Channel	Unit
PAP1-C2-20-35-OL	501034	Micro-Tip, CWT	20 mm	.035"	2.8 mm	5
PAP1-C2-25-35-OL	501036	Micro-Tip, CWT	25 mm	.035"	2.8 mm	5
PAP1-C2-30-35-OL	501039	Micro-Tip, CWT	30 mm	.035"	2.8 mm	5
PAP1-C5-20-35-OL	501046	Micro-Tip, monofil	20 mm	.035"	2.8 mm	5
PAP1-C5-30-35-OL	501047	Micro-Tip, monofil	30 mm	.035"	2.8 mm	5
PAP1-D2-30-21-OL	501853	Micro-Tip, CWT, Balloon Enteroscopy	30 mm	.021"	2.0 mm	2
PAP1-A2-20-35-OL	501028	Precut, CWT	20 mm	.035"	2.8 mm	5




PAP1-K2-06-35-OL	501060	Needle	6 mm	.035"	2.8 mm	5
------------------	--------	--------	------	-------	--------	---



PAP1-M2-20-35-OL	501840	B II, CWT	20 mm	.035"	2.8 mm	2
PAP1-N2-20-21-OL	501902	B II, CWT, Balloon Enteroscopy	20 mm	.021"	2.0 mm	2






Anatomical 3-D orientation
Based on systematic 11 o'clock orientation
the AXS_tome sphincterotomes also
follow the left frontal anatomy of
the distal choledochus.



AXS_tome Sphincterotomes for **easyCure**® (short and long wire | 3 lumen)



	Product Code	Order No.	Characteristics	Cutting Length	Accepted Guide Wire	ø Working Channel	Unit
	PAP1-C7-20-E35-OL	501930	Micro-Tip, CWT	20 mm	.035" ≥ 200 cm	3.8 mm	5
	PAP1-C7-30-E35-OL	501929	Micro-Tip, CWT	30 mm	.035" ≥ 200 cm	3.8 mm	5
	PAP1-C8-20-E35-OL	501931	Micro-Tip, monofil	20 mm	.035" ≥ 200 cm	3.8 mm	5
	PAP1-C8-30-E35-OL	501919	Micro-Tip, monofil	30 mm	.035" ≥ 200 cm	3.8 mm	5
	PAP1-F8-30-E35-OL	502017	Micro-Tip, monofil, insulated cutting wire (ProCut)	15 15 mm proximal insulated	.035" ≥ 200 cm	3.8 mm	5
	PAP1-K8-06-E35-OL	501920	Needle	6 mm	.035" ≥ 200 cm	3.8 mm	5
	PAP1-X8-30-35-260-OL	502286	Micro-Tip, monofil, preloaded gSlider	30 mm	.035" 260 cm	3.8 mm	5



easyCure®-Guide wire fixation



	Product Code	Order No.	Characteristics	Unit
	SPE1-Y7	501925	Guide wire fixation (instrument fork), for Olympus and Fujinon (from G5 series on)	5


GPS Wire Guide wires with hydrophilic tip



	Product Code	Order No.	Characteristics	Diameter	Length	Unit
	WIR1-S6-35-260	502517	Straight, Nitinol, insulated, green/white	.035"	260 cm	5
	WIR1-S6-35-450	502515	Straight, Nitinol, insulated, green/white	.035"	450 cm	5
	WIR1-J6-35-260	502518	J-tip, Nitinol, insulated, green/white	.035"	260 cm	5
	WIR1-J6-35-450	502516	J-tip, Nitinol, insulated, green/white	.035"	450 cm	5

ERCP-Cannulas




	Product Code	Order No.	Characteristics	Length	Accepted Guide Wire	ø Working Channel	Unit
	CAN1-A3-18-200-35	500337	Standard, 1 lumen, radiopaque marking	200 cm	.035"	2.0 mm	10
	CAN1-A5-18-200-35	500339	Standard, 1 lumen, metal tip	200 cm	.035"	2.8 mm	10
	CAN1-A6-18-200-35	502332	Standard, guide wire port, 1 lumen, metal tip, tapered (2.3 1.8 mm)	200 cm	.035"	2.8 mm	10
	CAN1-B3-18-200-35	500343	Filiform, 1 lumen, radiopaque marking	200 cm	.035"	2.0 mm	10
	CAN1-B3-18-260-35	501856	Filiform, 1 lumen, radiopaque marking, Balloon Enteroscopy	260 cm	.035"	2.0 mm	2
	CAN1-B4-18-200-35	502008	Filiform, 2 lumen, radiopaque marking, tapered (2.3 1.8 mm)	200 cm	.035"	2.8 mm	10
	CAN1-B4-18-200-E35	502007	Filiform, 2 lumen, radiopaque marking, for short and long wire technique, tapered (2.3 1.8 mm)	200 cm	.035" ≥ 200 cm	3.2 mm	2
	CAN1-D5-18-200-35	500348	Bottle neck, 1 lumen, metal tip	200 cm	.035"	2.0 mm	10


twist'n'CATCH-Stone extraction baskets




Product Code	Order No.	Characteristics	ø Basket	Length	Accepted Guide Wire	ø Working Channel	Unit
BAS1-K2-25-23-200	500027	Spherical, 130° rotation	35 mm	200 cm	-	2.8 mm	5
BAS1-K3-35-28-200	502285	Spherical, 130° rotation, for short and long wire technique	35 mm	200 cm	.035" ≥ 200 cm	3.2 mm	5



ATTENTION: medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.



For guide wires $\geq 200^*$



Maximum expansion power

*Figure: BAS1-K3-35-28-200

Dr. med. Martin Fährdrich Chief physician of interventional endoscopy at Dortmund Hospital



ROCKSTAR is easy to insert on account of its good flexibility in the distal area. Lithotripsy can be performed immediately without changing systems. That means material costs are kept within budget and the procedure is quicker. Overall, ROCKSTAR saves material and time so it is less expensive. In extreme situations, when performing lithotripsy on very hard stones, I recommend using a system with a metal coil.




The right tool for every situation

No matter how large the ‚rock‘ may be, Rockstar is fearless: with the new stone remover for lithotripsy you are ideally equipped for retrieving gallstones.

The guide wire lumen enables you to reach the location of the stone reliably. The innovative operating concept allows intuitive pulverisation of concretions if necessary.

ROCKSTAR Lithotriptor stone extraction basket




	Product Code	Order No.	Characteristics	ø Basket	Length	Accepted Guide Wire	ø Working Channel	Unit
	BAS3-K2-35-35-200	502486	twist'n'CATCH-Basket, integrated lithotripsy function	35 mm	200 cm	.035"	4.2 mm	2

ATTENTION: medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Stone extraction baskets



	Product Code	Order No.	Characteristics	ø Basket	Length	ø Working Channel	Unit
	BAS1-A2-20-23-200	500018	Hexagonal	20 mm	200 cm	2.8 mm	5
	BAS1-A2-25-23-200	500020	Hexagonal	25 mm	200 cm	2.8 mm	5
	BAS1-A2-30-23-200	500022	Hexagonal	30 mm	200 cm	2.8 mm	5
	BAS1-A2-35-23-200	500024	Hexagonal	35 mm	200 cm	2.8 mm	5
	BAS1-C1-20-18-260	500025	Half twisted, Balloon Enteroscopy	20 mm	260 cm	2.0 mm	2
	BAS1-C2-25-23-260	500026	Half twisted, Balloon Enteroscopy	25 mm	260 cm	2.8 mm	2




Not suitable for lithotripsy

ATTENTION: medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Stone extraction baskets (Reusable)

134°C

	Product Code	Order No.	Characteristics	ø Basket	Length	ø Working Channel	Unit
	twist'n'CATCH BAS2-K2-25-23-200	500120	Spherical, 130° rotation	35 mm	200 cm	2.8 mm	1
	BAS2-A2-25-23-200	500083	Hexagonal	25 mm	200 cm	2.8 mm	1
	BAS2-A2-30-23-200	500084	Hexagonal	30 mm	200 cm	2.8 mm	1
	BAS2-C1-15-18-200	500096	Half twisted	15 mm	200 cm	2.0 mm	1



Not suitable for lithotripsy

ATTENTION: medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

PD Dr. med. Alexander Dechêne Chief senior physician at the Gastroenterology and Hepatology Clinic, Essen University Hospital



The handling of a lithotripter should be uncomplicated because its use is relatively rare. rotaCRUSH is an easy-to-use one-piece system where operation is largely intuitive. Generally speaking, for lithotripsy systems it is absolutely essential to hold regular training in order to provide the endoscopy team with the necessary reliability in use.




The lithotripter for tough cases

Ready to use immediately without assembly. With the proven twist´n´CATCH basket you capture stones quickly and reliably. The crown head design of the metal coil ensures precision

force transmission to the stone and protects against basket detachment. The handle offers an integrated, additional Power² function for really tough cases.

rotaCRUSH Lithotripter





	Product Code	Order No.	Characteristics	ø Basket	Length	ø Working Channel	Unit
	LIT3-K2-30-35-200	502165	twist´n´CATCH-Basket, integrated catheter, spiral spring, handle	35 mm	200 cm	3.8 mm	1

ATTENTION: medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Lithotripter baskets





	Product Code	Order No.	Characteristics	ø Basket	Length	ø Working Channel	Unit
	twist'n'CATCH LIT1-K2-25-23-200	500826	Spherical, 130° rotation	35 mm	200 400 cm	2.8 mm	2
	LIT1-C3-20-23-200	500822	Half twisted	20 mm	200 400 cm	2.8 mm	2
	LIT1-C3-25-23-200	500823	Half twisted	25 mm	200 400 cm	2.8 mm	2
	LIT1-C3-30-23-200	500824	Half twisted	30 mm	200 400 cm	2.8 mm	2
	LIT1-C3-35-23-200	500825	Half twisted	35 mm	200 400 cm	2.8 mm	2

ATTENTION: medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Lithotripter baskets (Reusable)


134°C

	Product Code	Order No.	Characteristics	ø Basket	Length	ø Working Channel	Unit
	twist'n'CATCH LIT2-K2-25-23-200	500884	Spherical, 130° rotation	35 mm	200 400 cm	2.8 mm	1
	LIT2-C1-15-18-200	500871	Half twisted	15 mm	200 400 cm	2.0 mm	1
	LIT2-C3-20-23-200	500877	Half twisted	20 mm	200 400 cm	2.8 mm	1
	LIT2-C3-25-23-200	500878	Half twisted	25 mm	200 400 cm	2.8 mm	1
	LIT2-C3-30-23-200	500879	Half twisted	30 mm	200 400 cm	2.8 mm	1
	LIT2-C3-35-23-200	500880	Half twisted	35 mm	200 400 cm	2.8 mm	1

ATTENTION: medwork products of the BAS series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

KENA-Biliary stents (2 Flaps, centrally bent)




	Product Code	Order No.	Color	Diameter	Length between flaps	Total Length	Unit
	PR01-C2-07-05	502362	Grey	7 Fr	5 cm	6.5 cm	1
	PR01-C2-07-07	502368	Grey	7 Fr	7 cm	8.5 cm	1
	PR01-C2-07-09	502366	Grey	7 Fr	9 cm	10.5 cm	1
	PR01-C2-07-11	502367	Grey	7 Fr	11 cm	12.5 cm	1
	PR01-C2-07-13	502369	Grey	7 Fr	13 cm	14.5 cm	1
	PR01-C2-07-15	502370	Grey	7 Fr	15 cm	16.5 cm	1
	PR01-C2-07-17	502371	Grey	7 Fr	17 cm	18.5 cm	1
	PR01-C2-07-19	502372	Grey	7 Fr	19 cm	20.5 cm	1
	PR01-C2-08-05	502373	Green	8.5 Fr	5 cm	6.5 cm	1
	PR01-C2-08-07	502374	Green	8.5 Fr	7 cm	8.5 cm	1
	PR01-C2-08-09	502375	Green	8.5 Fr	9 cm	10.5 cm	1
	PR01-C2-08-11	502376	Green	8.5 Fr	11 cm	12.5 cm	1
	PR01-C2-08-13	502377	Green	8.5 Fr	13 cm	14.5 cm	1
	PR01-C2-08-15	502378	Green	8.5 Fr	15 cm	16.5 cm	1
	PR01-C2-08-17	502379	Green	8.5 Fr	17 cm	18.5 cm	1
	PR01-C2-08-19	502380	Green	8.5 Fr	19 cm	20.5 cm	1
	PR01-C2-10-05	502381	Blue	10 Fr	5 cm	6.5 cm	1
	PR01-C2-10-06	502419	Blue	10 Fr	6 cm	7.5 cm	1
	PR01-C2-10-07	502382	Blue	10 Fr	7 cm	8.5 cm	1
	PR01-C2-10-08	502420	Blue	10 Fr	8 cm	9.5 cm	1
	PR01-C2-10-09	502383	Blue	10 Fr	9 cm	10.5 cm	1
	PR01-C2-10-10	502421	Blue	10 Fr	10 cm	11.5 cm	1
	PR01-C2-10-11	502384	Blue	10 Fr	11 cm	12.5 cm	1
	PR01-C2-10-12	502422	Blue	10 Fr	12 cm	13.5 cm	1
	PR01-C2-10-13	502385	Blue	10 Fr	13 cm	14.5 cm	1
	PR01-C2-10-14	502423	Blue	10 Fr	14 cm	15.5 cm	1
	PR01-C2-10-15	502386	Blue	10 Fr	15 cm	16.5 cm	1
	PR01-C2-10-17	502387	Blue	10 Fr	17 cm	18.5 cm	1
	PR01-C2-10-19	502388	Blue	10 Fr	19 cm	20.5 cm	1

KENA-Biliary stents (2 Flaps, duodenal bent)




Product Code	Order No.	Color	Diameter	Length between flaps	Total Length	Unit
PR01-T2-07-05	502392	Grey	7 Fr	5 cm	6.5 cm	1
PR01-T2-07-07	502395	Grey	7 Fr	7 cm	8.5 cm	1
PR01-T2-07-09	502393	Grey	7 Fr	9 cm	10.5 cm	1
PR01-T2-07-11	502394	Grey	7 Fr	11 cm	12.5 cm	1
PR01-T2-07-13	502396	Grey	7 Fr	13 cm	14.5 cm	1
PR01-T2-07-15	502397	Grey	7 Fr	15 cm	16.5 cm	1
PR01-T2-07-17	502398	Grey	7 Fr	17 cm	18.5 cm	1
PR01-T2-07-19	502399	Grey	7 Fr	19 cm	20.5 cm	1
PR01-T2-08-05	502400	Green	8.5 Fr	5 cm	6.5 cm	1
PR01-T2-08-07	502401	Green	8.5 Fr	7 cm	8.5 cm	1
PR01-T2-08-09	502402	Green	8.5 Fr	9 cm	10.5 cm	1
PR01-T2-08-11	502403	Green	8.5 Fr	11 cm	12.5 cm	1
PR01-T2-08-13	502404	Green	8.5 Fr	13 cm	14.5 cm	1
PR01-T2-08-15	502405	Green	8.5 Fr	15 cm	16.5 cm	1
PR01-T2-08-17	502406	Green	8.5 Fr	17 cm	18.5 cm	1
PR01-T2-08-19	502407	Green	8.5 Fr	19 cm	20.5 cm	1
PR01-T2-10-05	502408	Blue	10 Fr	5 cm	6.5 cm	1
PR01-T2-10-06	502424	Blue	10 Fr	6 cm	7.5 cm	1
PR01-T2-10-07	502409	Blue	10 Fr	7 cm	8.5 cm	1
PR01-T2-10-08	502425	Blue	10 Fr	8 cm	9.5 cm	1
PR01-T2-10-09	502410	Blue	10 Fr	9 cm	10.5 cm	1
PR01-T2-10-10	502426	Blue	10 Fr	10 cm	11.5 cm	1
PR01-T2-10-11	502411	Blue	10 Fr	11 cm	12.5 cm	1
PR01-T2-10-12	502427	Blue	10 Fr	12 cm	13.5 cm	1
PR01-T2-10-13	502412	Blue	10 Fr	13 cm	14.5 cm	1
PR01-T2-10-14	502428	Blue	10 Fr	14 cm	15.5 cm	1
PR01-T2-10-15	502413	Blue	10 Fr	15 cm	16.5 cm	1
PR01-T2-10-17	502414	Blue	10 Fr	17 cm	18.5 cm	1
PR01-T2-10-19	502415	Blue	10 Fr	19 cm	20.5 cm	1

Preloaded **KENA**-Biliary stents (2 Flaps, centrally bent)



Product Code	Order No.	Color	Diameter	Length between flaps	Total Length	ø Working Channel	Unit
PR03-C2-08-05	502435	Green	8.5 Fr	5 cm	6.5 cm	4.2 mm	1
PR03-C2-08-07	502438	Green	8.5 Fr	7 cm	8.5 cm	4.2 mm	1
PR03-C2-08-09	502439	Green	8.5 Fr	9 cm	10.5 cm	4.2 mm	1
PR03-C2-08-11	502440	Green	8.5 Fr	11 cm	12.5 cm	4.2 mm	1
PR03-C2-08-13	502441	Green	8.5 Fr	13 cm	14.5 cm	4.2 mm	1
PR03-C2-08-15	502442	Green	8.5 Fr	15 cm	16.5 cm	4.2 mm	1
PR03-C2-08-17	502443	Green	8.5 Fr	17 cm	18.5 cm	4.2 mm	1
PR03-C2-08-19	502502	Green	8.5 Fr	19 cm	20.5 cm	4.2 mm	1
PR03-C2-10-05	502444	Blue	10 Fr	5 cm	6.5 cm	4.2 mm	1
PR03-C2-10-06	502445	Blue	10 Fr	6 cm	7.5 cm	4.2 mm	1
PR03-C2-10-07	502446	Blue	10 Fr	7 cm	8.5 cm	4.2 mm	1
PR03-C2-10-08	502447	Blue	10 Fr	8 cm	9.5 cm	4.2 mm	1
PR03-C2-10-09	502448	Blue	10 Fr	9 cm	10.5 cm	4.2 mm	1
PR03-C2-10-10	502449	Blue	10 Fr	10 cm	11.5 cm	4.2 mm	1
PR03-C2-10-11	502450	Blue	10 Fr	11 cm	12.5 cm	4.2 mm	1
PR03-C2-10-12	502451	Blue	10 Fr	12 cm	13.5 cm	4.2 mm	1
PR03-C2-10-13	502452	Blue	10 Fr	13 cm	14.5 cm	4.2 mm	1
PR03-C2-10-14	502453	Blue	10 Fr	14 cm	15.5 cm	4.2 mm	1
PR03-C2-10-15	502454	Blue	10 Fr	15 cm	16.5 cm	4.2 mm	1
PR03-C2-10-17	502455	Blue	10 Fr	17 cm	18.5 cm	4.2 mm	1
PR03-C2-10-19	502503	Blue	10 Fr	19 cm	20.5 cm	4.2 mm	1
PR03-C2-10-21	502504	Blue	10 Fr	21 cm	22.5 cm	4.2 mm	1
PR03-C2-10-23	502505	Blue	10 Fr	23 cm	24.5 cm	4.2 mm	1
PR03-C2-10-25	502506	Blue	10 Fr	25 cm	26.5 cm	4.2 mm	1


Preloaded **KENA**-Biliary stents (2 Flaps, duodenal bent)



Product Code	Order No.	Color	Diameter	Length between flaps	Total Length	ø Working Channel	Unit
PR03-T2-08-05	502436	Green	8.5 Fr	5 cm	6.5 cm	4.2 mm	1
PR03-T2-08-07	502456	Green	8.5 Fr	7 cm	8.5 cm	4.2 mm	1
PR03-T2-08-09	502457	Green	8.5 Fr	9 cm	10.5 cm	4.2 mm	1
PR03-T2-08-11	502458	Green	8.5 Fr	11 cm	12.5 cm	4.2 mm	1
PR03-T2-08-13	502459	Green	8.5 Fr	13 cm	14.5 cm	4.2 mm	1
PR03-T2-08-15	502460	Green	8.5 Fr	15 cm	16.5 cm	4.2 mm	1
PR03-T2-08-17	502461	Green	8.5 Fr	17 cm	18.5 cm	4.2 mm	1
PR03-T2-08-19	502507	Green	8.5 Fr	19 cm	20.5 cm	4.2 mm	1
PR03-T2-10-05	502462	Blue	10 Fr	5 cm	6.5 cm	4.2 mm	1
PR03-T2-10-06	502463	Blue	10 Fr	6 cm	7.5 cm	4.2 mm	1
PR03-T2-10-07	502464	Blue	10 Fr	7 cm	8.5 cm	4.2 mm	1
PR03-T2-10-08	502465	Blue	10 Fr	8 cm	9.5 cm	4.2 mm	1
PR03-T2-10-09	502466	Blue	10 Fr	9 cm	10.5 cm	4.2 mm	1
PR03-T2-10-10	502467	Blue	10 Fr	10 cm	11.5 cm	4.2 mm	1
PR03-T2-10-11	502468	Blue	10 Fr	11 cm	12.5 cm	4.2 mm	1
PR03-T2-10-12	502469	Blue	10 Fr	12 cm	13.5 cm	4.2 mm	1
PR03-T2-10-13	502470	Blue	10 Fr	13 cm	14.5 cm	4.2 mm	1
PR03-T2-10-14	502471	Blue	10 Fr	14 cm	15.5 cm	4.2 mm	1
PR03-T2-10-15	502472	Blue	10 Fr	15 cm	16.5 cm	4.2 mm	1
PR03-T2-10-17	502473	Blue	10 Fr	17 cm	18.5 cm	4.2 mm	1
PR03-T2-10-19	502508	Blue	10 Fr	19 cm	20.5 cm	4.2 mm	1
PR03-T2-10-21	502509	Blue	10 Fr	21 cm	22.5 cm	4.2 mm	1
PR03-T2-10-23	502510	Blue	10 Fr	23 cm	24.5 cm	4.2 mm	1
PR03-T2-10-25	502511	Blue	10 Fr	25 cm	26.5 cm	4.2 mm	1




Tips tapered on both sides

The unique bilateral tip shaping of the duodenal Kena prostheses facilitates placement and minimises the perforation risk. Alternatively, the prostheses can be laid across a bifurcation with the bend first and thus be anchored in place. This ensures a better hold and reduces the risk of migration.




KENA-Guide catheters, -Pusher and -Insertion kit



	Product Code	Order No.	Characteristics	Diameter	Length	Accepted Guide Wire	ø Working Channel	Unit
	PR01-K1-05-220-35	502491	Guide catheter, grey, for 8.5 and 10 Fr Biliary stents	5 Fr	220 cm	.035"	2.0 mm	2
	PR01-P2-07-180	502494	Pusher, white, for 7 Fr Biliary stents	7 Fr	180 cm	.035"	2.8 mm	2
	PR01-P1-09-180	502492	Pusher, white, for 8.5 and 10 Fr Biliary stents	9 Fr	180 cm	.035"	3.2 mm	2
	PR01-S1-05-220-35	502493	Insertion kit for 8.5 and 10 Fr Biliary stents, LL-connection, distal radiopaque marking	9 Fr	220 cm	.035"	4.2 mm	2


ERCP

Stone extraction balloons



Product Code	Order No.	Characteristics	ø Balloon	Length	Accepted Guide Wire	ø Working Channel	Unit
BAL1-B2-16-75-200-35	500005	Distal radiopaque marking, 2 lumen, tapered (7 5 Fr)	16 mm	200 cm	.035"	2.8 mm	5
BAL1-C2-16-07-200-35	500010	Distal radiopaque marking, 2 lumen	16 mm	200 cm	.035"	2.8 mm	5
BAL1-C3-18-07-200-35	500011	Distal radiopaque marking, 3 lumen	16 mm	200 cm	.035"	2.8 mm	5
BAL1-C3-18-07-200E35	501916	Distal radiopaque marking, tapered, 3 lumen, for short and long wire technique	16 mm	200 cm	.035" ≥ 200 cm	3.2 mm	2

ERCP Kit



Product Code	Order No.	Characteristics	Unit
CAN3-STANDARD-1	500427	Disposable ERCP kit including: 1 Pair disposable gloves, medium, powder- and latex-free 1 Cup, 210 ml, white transparent 1 Cup, 210 ml, blue transparent 10 Compresses, 10 x 10 cm, 8-fold 1 Disposable syringe, 10 ml, Luer-Lock 1 Disposable syringe, 20 ml, Luer-Lock 1 Disposable ERCP towel, 75 x 90 cm, with 4 transparent insertion bags (MIS1-TOWEL-4)	20

Dilation balloons


Working Length 230 cm



	Product Code	Order No.	Characteristics	ø Balloon	Length Balloon	Accepted Guide Wire	ø Working Channel	Unit
	DIL1-A1-06-40-23-35	502268	2 Radiopaque markings, tapered catheter (2,3 1,8 mm), for short and long wire technique	6 mm	40 mm	.035" ≥ 200 cm	2.8 mm	2
	DIL1-A1-08-40-23-35	502269	2 Radiopaque markings, tapered catheter (2,3 1,8 mm), for short and long wire technique	8 mm	40 mm	.035" ≥ 200 cm	2.8 mm	2
	DIL1-A1-10-40-23-35	502270	2 Radiopaque markings, tapered catheter (2,3 1,8 mm), for short and long wire technique	10 mm	40 mm	.035" ≥ 200 cm	2.8 mm	2
	DIL1-A1-12-60-23-35	502271	2 Radiopaque markings, for short and long wire technique	12 mm	60 mm	.035" ≥ 200 cm	2.8 mm	2
	DIL1-A1-14-60-23-35	502272	2 Radiopaque markings, for short and long wire technique	14 mm	60 mm	.035" ≥ 200 cm	2.8 mm	2
	DIL1-A1-16-60-23-35	502273	2 Radiopaque markings, for short and long wire technique	16 mm	60 mm	.035" ≥ 200 cm	2.8 mm	2
	DIL1-A1-18-60-23-35	502274	2 Radiopaque markings, for short and long wire technique	18 mm	60 mm	.035" ≥ 200 cm	2.8 mm	2
	DIL1-A1-20-60-23-35	502275	2 Radiopaque markings, for short and long wire technique	20 mm	60 mm	.035" ≥ 200 cm	3.2 mm	2


Insufflation device



	Product Code	Order No.	Characteristics	Unit
	DIL1-B1-60	502278	60 ml volume, 12 Atm maximum pressure	1


JAZZ Endoscope valve kit for FUJIFILM 700 series



	Product Code	Order No.	Characteristics	Unit
	VAL1-G7-50	502562	Disposable endoscope valve kit for Fujifilm endoscopes (G7 series), compatible with ED-580XT and EI-580BT, consisting of: air/water valve, suction valve and biopsy valve FROG	50

FROG Biopsy valves



	Product Code	Order No.	Characteristics	Unit
	VAL1-F1-100	502484	For Olympus and Fujifilm endoscopes	100
	VAL1-P1-100	502485	For Pentax endoscopes	100

GPS Wire Guide wires with hydrophilic tip



	Product Code	Order No.	Characteristics	Diameter	Length	Unit
	WIR1-S6-35-260	502517	Straight, Nitinol, insulated, green/white	.035"	260 cm	5
	WIR1-S6-35-450	502515	Straight, Nitinol, insulated, green/white	.035"	450 cm	5
	WIR1-J6-35-260	502518	J-tip, Nitinol, insulated, green/white	.035"	260 cm	5
	WIR1-J6-35-450	502516	J-tip, Nitinol, insulated, green/white	.035"	450 cm	5

Guide wires



	Product Code	Order No.	Characteristics	Diameter	Length	Unit
	WIR1-E1-72-260	501780	Straight, steel, flexible tip, Eder-Püstow	.072"	260 cm	2
	WIR1-F1-18-400	501781	Straight, Nitinol, coated	.018"	400 cm	2
	WIR1-F1-21-400	501782	Straight, Nitinol, coated	.021"	400 cm	2
	WIR1-F1-21-500	501911	Straight, Nitinol, coated	.021"	500 cm	2
	WIR1-F5-25-400	501789	Straight, Nitinol, insulated, black/white	.025"	400 cm	2
	WIR1-F5-35-400	501792	Straight, Nitinol, insulated, black/white	.035"	400 cm	2
	WIR1-F7-35-260	501837	Straight, Nitinol, coated, extra stiff	.035"	260 cm	2
	WIR1-F7-35-400	501838	Straight, Nitinol, coated, extra stiff	.035"	400 cm	2

1-9

Lithotripter baskets



	Product Code	Order No.	Characteristics	ø Basket	Length	ø Working Channel	Unit
	twist'n'CATCH	500826	Spherical, 130° rotation	35 mm	200 400 cm	2,8 mm	2
	LIT1-K2-25-23-200						
	LIT1-C3-20-23-200	500822	Half twisted	20 mm	200 400 cm	2,8 mm	2
	LIT1-C3-25-23-200	500823	Half twisted	25 mm	200 400 cm	2,8 mm	2
	LIT1-C3-30-23-200	500824	Half twisted	30 mm	200 400 cm	2,8 mm	2
	LIT1-C3-35-23-200	500825	Half twisted	35 mm	200 400 cm	2,8 mm	2

CAUTION: medwork products of the LIT series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Lithotripter springs



	Product Code	Order No.	Characteristics	Length	ø Working Channel	Unit
	LIT1-G4-26-185	502304	Spiral spring, transendoscopic, bile duct	185 cm	2,8 mm	2
	LIT1-N4-38-090	502305	Spiral spring, emergency, bile duct	90 cm	-	2

CAUTION: medwork products of the LIT series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Handle (Reusable)

134°C


	Product Code	Order No.	Characteristics	Unit
	LIT2-D-HAN	500881	Handle, pull-through technique	1

1-15

CAUTION: medwork products of the LIT series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.

Preloaded **KENA** Biliary stents (2 flaps, centre bend) **easycure**

For short wire technique



Product Code	Order No.	Color	Diameter	Length Between Flaps	Total Length	ø Working Channel	Unit
PR03-C3-08-05	502652	Green	8.5 Fr	5 cm	6,5 cm	4,2 mm	1
PR03-C3-08-07	502653	Green	8.5 Fr	7 cm	8,5 cm	4,2 mm	1
PR03-C3-08-09	502654	Green	8.5 Fr	9 cm	10,5 cm	4,2 mm	1
PR03-C3-08-11	502655	Green	8.5 Fr	11 cm	12,5 cm	4,2 mm	1
PR03-C3-08-13	502656	Green	8.5 Fr	13 cm	14,5 cm	4,2 mm	1
PR03-C3-08-15	502657	Green	8.5 Fr	15 cm	16,5 cm	4,2 mm	1
PR03-C3-08-17	502658	Green	8.5 Fr	17 cm	18,5 cm	4,2 mm	1
PR03-C3-08-19	502659	Green	8.5 Fr	19 cm	20,5 cm	4,2 mm	1
PR03-C3-10-05	502660	Blue	10 Fr	5 cm	6,5 cm	4,2 mm	1
PR03-C3-10-06	502661	Blue	10 Fr	6 cm	7,5 cm	4,2 mm	1
PR03-C3-10-07	502662	Blue	10 Fr	7 cm	8,5 cm	4,2 mm	1
PR03-C3-10-08	502663	Blue	10 Fr	8 cm	9,5 cm	4,2 mm	1
PR03-C3-10-09	502664	Blue	10 Fr	9 cm	10,5 cm	4,2 mm	1
PR03-C3-10-10	502665	Blue	10 Fr	10 cm	11,5 cm	4,2 mm	1
PR03-C3-10-11	502666	Blue	10 Fr	11 cm	12,5 cm	4,2 mm	1
PR03-C3-10-12	502667	Blue	10 Fr	12 cm	13,5 cm	4,2 mm	1
PR03-C3-10-13	502668	Blue	10 Fr	13 cm	14,5 cm	4,2 mm	1
PR03-C3-10-14	502669	Blue	10 Fr	14 cm	15,5 cm	4,2 mm	1
PR03-C3-10-15	502670	Blue	10 Fr	15 cm	16,5 cm	4,2 mm	1
PR03-C3-10-17	502671	Blue	10 Fr	17 cm	18,5 cm	4,2 mm	1
PR03-C3-10-19	502672	Blue	10 Fr	19 cm	20,5 cm	4,2 mm	1
PR03-C3-10-21	502673	Blue	10 Fr	21 cm	22,5 cm	4,2 mm	1
PR03-C3-10-23	502674	Blue	10 Fr	23 cm	24,5 cm	4,2 mm	1
PR03-C3-10-25	502675	Blue	10 Fr	25 cm	26,5 cm	4,2 mm	1

Preloaded **KENA** Biliary stents (2 flaps, proximal bend) **easyCure®**

For short wire technique



	Product Code	Order No.	Color	Diameter	Length Between Flaps	Total Length	ø Working Channel	Unit
	PR03-T3-08-05	502676	Green	8.5 Fr	5 cm	6,5 cm	4,2 mm	1
	PR03-T3-08-07	502677	Green	8.5 Fr	7 cm	8,5 cm	4,2 mm	1
	PR03-T3-08-09	502678	Green	8.5 Fr	9 cm	10,5 cm	4,2 mm	1
	PR03-T3-08-11	502679	Green	8.5 Fr	11 cm	12,5 cm	4,2 mm	1
	PR03-T3-08-13	502680	Green	8.5 Fr	13 cm	14,5 cm	4,2 mm	1
	PR03-T3-08-15	502681	Green	8.5 Fr	15 cm	16,5 cm	4,2 mm	1
	PR03-T3-08-17	502682	Green	8.5 Fr	17 cm	18,5 cm	4,2 mm	1
	PR03-T3-08-19	502683	Green	8.5 Fr	19 cm	20,5 cm	4,2 mm	1
	PR03-T3-10-05	502684	Blue	10 Fr	5 cm	6,5 cm	4,2 mm	1
	PR03-T3-10-06	502685	Blue	10 Fr	6 cm	7,5 cm	4,2 mm	1
	PR03-T3-10-07	502686	Blue	10 Fr	7 cm	8,5 cm	4,2 mm	1
	PR03-T3-10-08	502687	Blue	10 Fr	8 cm	9,5 cm	4,2 mm	1
	PR03-T3-10-09	502688	Blue	10 Fr	9 cm	10,5 cm	4,2 mm	1
	PR03-T3-10-10	502689	Blue	10 Fr	10 cm	11,5 cm	4,2 mm	1
	PR03-T3-10-11	502690	Blue	10 Fr	11 cm	12,5 cm	4,2 mm	1
	PR03-T3-10-12	502691	Blue	10 Fr	12 cm	13,5 cm	4,2 mm	1
	PR03-T3-10-13	502692	Blue	10 Fr	13 cm	14,5 cm	4,2 mm	1
	PR03-T3-10-14	502693	Blue	10 Fr	14 cm	15,5 cm	4,2 mm	1
	PR03-T3-10-15	502694	Blue	10 Fr	15 cm	16,5 cm	4,2 mm	1
	PR03-T3-10-17	502695	Blue	10 Fr	17 cm	18,5 cm	4,2 mm	1
	PR03-T3-10-19	502696	Blue	10 Fr	19 cm	20,5 cm	4,2 mm	1
	PR03-T3-10-21	502697	Blue	10 Fr	21 cm	22,5 cm	4,2 mm	1
	PR03-T3-10-23	502698	Blue	10 Fr	23 cm	24,5 cm	4,2 mm	1
	PR03-T3-10-25	502699	Blue	10 Fr	25 cm	26,5 cm	4,2 mm	1

ESD ACCESSORIES

Model	Outer diameter (mm)	Inner diameter of tip (mm)	Tip length (mm)	Accessory guide	Applicable endoscope
-------	------------------------	-------------------------------	--------------------	-----------------	----------------------

ST hood for ESD and POEM

DH-15GR	12.2	7.0	8.3	Yes	EG-760Z, EC-740T
DH-16CR	16.0	7.0	8.3	Yes	EC-760R, EC-760ZP
DH-28GR	11.8	8.0	7.0	Without	EG-760R, EG-760Z, EC-740T
DH-29CR	13.0	8.0	7.0	Without	EG-760CT, EC-760P
DH-30CR	14.8	8.0	7.0	Without	EC-760R, EC-760ZP
DH-33GR	NEW 11.8	7.0	7.0	Yes	EG-760Z, EC-740T
DH-34CR	NEW 13.6	7.0	7.0	Yes	EC-760R, EC-760ZP

Observation hood

DH-11GZ	13.2		2.5		EG-590ZW
DH-13CZ	15.7		3.0		EC-760ZP-VM/VL, EC-760R-VM/VI/VL, EC-600WM, WI, WL, EC-600ZW M/L, EC-590WM4, WI4, WL4, EC-590ZW3 M/L, EC-530WM3, WI3, WL3
DH-14EN	10.8		1.5		EG-760Z, EG-600ZW
DH-17EN	11.5		1.5		EI-580BT, EN-580T



Eleview™

Submucosal injection agent designed for easy and safe resection procedures.

Manufactured by Cosmo Technologies Ltd., Ireland. Not for sale in Canada.

Model	Applicable endoscope	Remarks	THIRD PARTY PRODUCTS
-------	----------------------	---------	----------------------

Top Endoscopic Hoods


SHS	EG-760R	Soft & transparent type
SHM	EG-760Z, EG-760CT, EC-740T	Soft & transparent type
SHL	EC-760R, EC-760ZP, EC-760P	Soft & transparent type
SHM-B	EG-760Z, EG-760CT, EC-740T	Soft & black type
SHL-B	EC-760R, EC-760ZP, EC-760P	Soft & black type
SHLL	EC-760Z	Soft & transparent type

Packaging unit: 5 pcs/box

Manufactured by Top Corporation, Japan. Not for sale in Canada.

Model	Product Description	THIRD PARTY PRODUCTS
-------	---------------------	----------------------

Monitors

Radiance 27" ULTRA 	 Image: Radiance 27"	27" HD type LCD monitor with Ultra bright LED Backlight ▶ High-Definition, Color Correction Technology (CCT), Full Multi-Modality, Gorilla Glass front panel ▶ 8.9 kg, 678 mm (W) x 445 mm (H) x 84 mm (D) ▶ Input signal: HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA ▶ Output signal: HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite
Radiance 26" G2 HB LED 		26" HD type High Brightness LCD monitor with LED Backlight ▶ High-Definition, Color Correction Technology (CCT), Full Multi-Modality ▶ 8.2 kg, 673 mm (W) x 418 mm (H) x 88 mm (D) ▶ Input signal: HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA ▶ Output signal: HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite
Radiance 24" G2 LED 		24" HD type LCD monitor with LED Backlight ▶ High-Definition, Color Correction Technology (CCT), Full Multi-Modality ▶ 7.1 kg, 597 mm (W) x 401 mm (H) x 100 mm (D) ▶ Input signal: HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA ▶ Output signal: HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite
Radiance 19" HD 		19" HD type LCD monitor ▶ High-Definition, Full Multi-Modality ▶ 6.8 kg, 465 mm (W) x 400 mm (H) x 98 mm (D) ▶ Input signal: HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, composite, VGA ▶ Output signal: HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite
Including FUJIFILM BIOS for the best performance		
EndoVue 19" 		19" HD type LCD monitor with LED Backlight ▶ 4.2 kg, 464.8 mm (W) x 396.2 mm (H) x 99 mm (D) ▶ Input Signal: DVI-D, HD-SDI, (HD-)RGBS, (HD-)YPbPr, VGA, S-Video, Composite

Manufactured by NDS Surgical Imaging, LLC, USA



5 ampoules/pouch
(10ml/ampoule)

Model	Applicable endoscope	Remarks	THIRD PARTY PRODUCTS
-------	----------------------	---------	----------------------

SHF-020	EG-760Z, EC-740T	Hard and transparent type
SHF-035	EG-760CT, EC-760P	Hard and transparent type
SHF-045	EC-760R	Hard and transparent type
SHF-050	EC-760ZP	Hard and transparent type
SHF-055	EC-760Z	Hard and transparent type

ESD ACCESSORIES

Model	Identifier	Knife length	Working length	Ø Suitable working channel	Remarks
		(mm)	(mm)		

FlushKnife

	DK2620J	- B15S -	1.5	2,000	2.8	Ball tip, Slim Type
		- B20S -	2.0			
		- B25S -	2.5			
		- B30S -	3.0			
	DK2623JS	- B15 -	1.5	2,300	2.8	Ball tip, Long slim type
		- B20 -	2.0			
		- N10S -	1.0			
		- N15S -	1.5			
	DK2620J	- N20S -	2.0	2,000	2.8	Needle tip, Slim Type
		- N25S -	2.5			
		- N30S -	3.0			
		- N15 -	1.5			
	DK2623JS	- N20 -	2.0	2,300	2.8	Needle tip, Long slim type

Model	Identifier	Knife length	Working length	Ø Suitable working channel	Remarks
		(mm)	(mm)		

ClutchCutter

	DP2618DT - 35 -	- 35 -	3.5	1,800	2.8
	DP2618DT - 50 -	- 50 -	5.0	1,800	2.8



MANTA Injection needles



	Product Code	Order No.	Characteristics	Needle ø Length	Length	ø Working Channel	Unit
	INJ1-A3-05-5-18-240	502488	Latching mechanism, distal needle guide, suitable with Histoacryl	0,5 mm (25 Gauge) 5 mm	240 cm	2,0 mm	5
	INJ1-A3-07-5-18-240	502487	Latching mechanism, distal needle guide, suitable with Histoacryl	0,7 mm (23 Gauge) 5 mm	240 cm	2,0 mm	5

For many injection needles, applications in the terminal ileum or caecum constitute borderline cases where numerous instruments can no longer be activated reliably. Even 90° angulation of the injection needle by the Albarran lever in a duodenoscope often leads to failure of the instrumentation. In such a case, Manta ensures reliable retraction and extension of the needle, even in sharply angulated device positions, and prevents feared damage to the working channel.



Polypectomy


Great effect on small polyps

Less effort and perfect cutting performance - that's what CrossSnare ZERO offers you thanks to its unique reinforced spiral tube design.


Equipped with extra thin cutting wire and available in two different sizes, CrossSnare ZERO is the perfect choice for polyps up to 9 mm and supports you optimally in implementing the ESGE guideline*.

CrossSnare ZERO Polypectomy snares for cold snaring



	Product Code	Order No.	Characteristics	ø Snare	Length	ø Working Channel	Unit
	POL1-Z2-07-23-220	502569	Kite, CWT	7 mm	220 cm	2,8 mm	10
	POL1-Z2-12-23-220	502570	Kite, CWT	12 mm	220 cm	2,8 mm	10

Polypectomy



More than just a cold snare - a hybrid comprising HOT and COLD

Over 90% of all polyps are smaller than 10 mm. Especially a number of small, flat lesions in the caecum or right colon constitute special challenges for the instrumentation. CrossSnare is the first hybrid snare that has been specially developed for the majority of polyps of clinical relevance in

terms of ablation. The large opening angle, even if the snare is partially retracted, ensures optimal snare adjustment to match the size of the polyps between 3 mm and 15 mm. Compacted Wire Technology provides tissue-sparing cutting properties - both with and without HF current.

CrossSnare Polypectomy snares for Cold and HF snaring

	Product Code	Order No.	Characteristics	ø Snare	Length	ø Working Channel	Unit
	POL1-X1-10-23-220-OL	502308	Kite, CWT	10 mm	220 cm	2,8 mm	10
	POL1-X1-15-23-220-OL	502309	Kite, CWT	15 mm	220 cm	2,8 mm	10

Polypectomy





Figure: POL1-C5-22-23-220-OL

ResectionMaster Polypectomy snares for EMR



	Product Code	Order No.	Characteristics	ø Snare	Length	ø Working Channel	Unit
	POL1-C5-10-23-220-OL	502293	Oval, flat metal band	10 mm	220 cm	2,8 mm	10
	POL1-C5-15-23-220-OL	502525	Oval, flat metal band	15 mm	220 cm	2,8 mm	10
	POL1-C5-22-23-220-OL	501896	Oval, flat metal band	22 mm	220 cm	2,8 mm	10
	POL1-C5-30-23-220-OL	501222	Oval, flat metal band	30 mm	220 cm	2,8 mm	10

Polypectomy





Figure: POL1-H4-15-27-230-OL

SnareInflator Polypectomy snares with integrated injection needle








Product Code	Order No.	Characteristics	ø Snare	Needle Length	Length	ø Working Channel	Unit
POL1-H4-15-27-230-OL	502346	Oval, CWT	15 mm	5 mm	230 cm	3,2 mm	5
POL1-H3-22-27-230-OL	501231	Oval, monofil	22 mm	5 mm	230 cm	3,2 mm	5
POL1-H3-30-27-230-OL	501948	Oval, monofil	30 mm	5 mm	230 cm	3,2 mm	5

	Product Code	Order No.	Characteristics	Needle ø Length	Length	ø Working Channel	Unit
	INJ1-A2-05-3-23-230	502597	Latching mechanism, distal needle guide	0.5 mm (25 Gauge) 3 mm	230 cm	2.8 mm	10
	INJ1-A2-07-3-23-230	502598	Latching mechanism, distal needle guide	0.7 mm (23 Gauge) 3 mm	230 cm	2.8 mm	10
	INJ1-A2-05-4-23-180	502333	Latching mechanism, distal needle guide	0.5 mm (25 Gauge) 4 mm	180 cm	2.8 mm	10
	INJ1-A2-05-4-23-230	502334	Latching mechanism, distal needle guide	0.5 mm (25 Gauge) 4 mm	230 cm	2.8 mm	10
	INJ1-A2-07-4-23-180	502335	Latching mechanism, distal needle guide	0.7 mm (23 Gauge) 4 mm	180 cm	2.8 mm	10
	INJ1-A2-07-4-23-230	502336	Latching mechanism, distal needle guide	0.7 mm (23 Gauge) 4 mm	230 cm	2.8 mm	10
	INJ1-A2-07-5-23-180	502337	Latching mechanism, distal needle guide	0.7 mm (23 Gauge) 5 mm	180 cm	2.8 mm	10
	INJ1-A2-07-5-23-230	502338	Latching mechanism, distal needle guide	0.7 mm (23 Gauge) 5 mm	230 cm	2.8 mm	10
	INJ1-A2-09-5-23-180	502433	Latching mechanism, distal needle guide	0.9 mm (21 Gauge) 5 mm	180 cm	2.8 mm	10
	INJ1-A2-09-5-23-230	502434	Latching mechanism, distal needle guide	0.9 mm (21 Gauge) 5 mm	230 cm	2.8 mm	10


	Product Code	Order No.	Characteristics	Needle ø length	Length	ø Working Channel	Unit
	INJ1-A1-05-5-18-180	500745	Latching mechanism, distal needle guide	0.5 mm (25 Gauge) 5 mm	180 cm	2.0 mm	10
	INJ1-A1-07-5-18-180	500752	Latching mechanism, distal needle guide	0.7 mm (23 Gauge) 5 mm	180 cm	2.0 mm	10
	INJ1-A1-07-8-23-230	500771	Latching mechanism, distal needle guide	0.7 mm (23 Gauge) 8 mm	180 cm	2.3 mm	10

OCTOPUS Foreign body retrieval net



	Product Code	Order No.	Characteristics	Volume Net	Length	ø Working Channel	Unit
	FOR1-N1-20-180	502477	Hexagonal, framed, paediatric	3 cm ³	180 cm	2,0 mm	5
	FOR1-N1-30-220	502479	Oval, framed	6 cm ³	220 cm	2,8 mm	5
	FOR1-N1-40-220	502480	Reservoir function, framed	8 cm ³	220 cm	2,8 mm	5
	FOR1-N2-30-220	502481	Kite, framed	6 cm ³	220 cm	2,8 mm	5
	FOR1-N3-30-220	502482	Kite, framed, rotatable	6 cm ³	220 cm	2,8 mm	5


Owing to a new reservoir design Octopus enables easy and reliable recovery of multiple objects simultaneously. The 1:1 rotatability facilitates manoeuvring during foreign body recovery. The extremely sturdy net reliably maintains its position and does not slip.



Polypectomy


Foreign body forceps



	Product Code	Order No.	Characteristics	Opening Width	ø Tube	Length	ø Working Channel	Unit
	FOR1-A1-23-230	502185	Alligator, metal spring, coated	9 mm	2,3 mm	230 cm	2,8 mm	5
	FOR1-A2-23-230	502284	Hybrid jaw configuration alligator/ rat tooth, metal spring, coated	9 mm	2,3 mm	230 cm	2,8 mm	10
	FOR1-A2-23-260	502548	Hybrid jaw configuration alligator/ rat tooth, metal spring, coated balloon enteroscopy	8 mm	2,3 mm	260 cm	2,8 mm	5
	FOR1-P1-18-230	502549	Pelican, coated, paediatric	9 mm	1,8 mm	230 cm	2,0 mm	5
	FOR1-P1-23-230	502547	Pelican, coated	15,5 mm	2,3 mm	230 cm	2,8 mm	5
	FOR1-H3-23-230	501946	Three-arm grasping forceps, PTFE catheter	-	2,3 mm	230 cm	2,8 mm	5


Protection cap



	Product Code	Order No.	Characteristics	Unit
	WS0420000*	500944	Foreign body protection cap	10

* Only available in Europe

KOALA Polyptrap




Product Code	Order No.	Characteristics	Unit
TRP1-S2	502499	2 integrated strainers, white	25
TRP1-S4	502536	2 integrated strainers, white and 2 additional strainers, grey	25

Minimised contamination risk

On account of the two integrated strainers of the new Koala it is no longer necessary to deposit a contaminated strainer outside the system. That minimises contamination with biomaterial, and hence the infection risk.



Polyptraps




Product Code	Order No.	Characteristics	Unit
MIS1-POLYPTRAP	500941	4 chambers, connecting tubes to suction pump endoscope	10
MIS1-POLYPTRAP50	502194	4 chambers, connecting tubes to suction pump endoscope	50

FROG Biopsy valves



Product Code	Order No.	Characteristics	Unit
VAL1-F1-100	502484	For Olympus and Fujifilm endoscopes	100
VAL1-P1-100	502485	For Pentax endoscopes	100

Spray catheters



Product Code	Order No.	Characteristics	Length	ø Working Channel	Unit
SPE1-B6-18-120	501847	Spray catheter	120 cm	2,0 mm	10
SPE1-C6-18-220	501718	Spray catheter	220 cm	2,0 mm	10

Cytology brushes




	Product Code	Order No.	Characteristics	ø Brush	Length	ø Working Channel	Unit
	WS0372512*	502159	Bronchoscopy, handle, 1 lumen, distal protection ball	2,5 mm	120 cm	2,0 mm	5
	WS0373016*	502161	Gastroscopy, handle, 1 lumen, distal protection ball	3,0 mm	160 cm	2,8 mm	5
	WS0373022*	502160	Colonoscopy, handle, 1 lumen, distal protection ball	3,0 mm	220 cm	2,8 mm	5

* Only available in Europe

Colonoscopy examination pants






	Product Code	Order No.	Characteristics	Unit
	MIS1-COLO-PANTS1	502276	Colonoscopy examination pants, dark blue, fleece, including elastic, slit, bermuda (width 61 cm length 65 cm)	100
	MIS1-COLO-PANTS	500920	Colonoscopy examination pants, dark blue, fleece, including elastic, slit with hook-and-loop cover, prolonged waistline, bermuda (width 70 cm length 60 cm)	100
	MIS1-COLO-PANTS-XXL	500921	Colonoscopy examination pants, dark blue, fleece, including elastic, slit with hook-and-loop cover, prolonged waistline, bermuda (width 85 cm length 70 cm)	10

Accessories

Miscellaneous



	Product Code	Order No.	Characteristics	Unit
	MIS1-COLO-PANTS1	502276	Colonoscopy examination pants, dark blue, fleece, including elastic, slit, bermuda (width 61 cm length 65 cm)	100
	MIS1-COLO-PANTS	500920	Colonoscopy examination pants, dark blue, fleece, including elastic, slit with hook-and-loop cover, prolonged waistline, bermuda (width 70 cm length 60 cm)	100
	MIS1-COLO-PANTS-XXL	500921	Colonoscopy examination pants, dark blue, fleece, including elastic, slit with hook-and-loop cover, prolonged waistline, bermuda (width 85 cm length 70 cm)	10
	MIS1-TOWEL-3*	500951	ERCP towel, sterile, with 3 transparent insertion bags (30 x 26 cm), 75 x 90 cm	10
	MIS1-TOWEL-4*	500952	ERCP towel, sterile, with 4 transparent insertion bags (30 x 20 cm), 75 x 90 cm	10

* Only available in Europe

Product clip (Reusable)


134°C

	Product Code	Order No.	Characteristics	Unit
	MIS2-PRO-CLIP	500975	Reusable product clip for fixation and reprocessing of flexible accessories, autoclavable up to 134° C, opening width 0-10 mm	20


Buried-Bumper-Therapy

FLAMINGO Kit




	Product Code	Order No.	Characteristics	Cutting Length	Length	Accepted Guide Wire	Unit
	POL3-P1-40-30-30	502490	Flamingo kit consisting of:				1
			1 Flamingo for 15 Fr PEG probe, 2 lumen, CWT, rotatable	6 10 mm	30 cm	.035"	
			1 Foreign body forceps, alligator/rat tooth	-	230 cm	-	
			1 Guide wire, hydrophilic tip	-	120 cm	-	


PUSH TECHNIQUE



1 Exploration with the PEG probe using a guide wire



2 Cutting the PEG probe free with the **Flamingo** by making a star-shaped incision








3 Removal of the PEG probe from the stomach with the foreign body forceps

Biopsy

EndoBite² Biopsy forceps



	Product Code	Order No.	Characteristics	ø Jaws	Length	ø Working Channel	Unit
	BI01-C4-18-120	502004	Oval jaws, coated, bronchoscopy	1,8 mm	120 cm	2,0 mm	10
	BI01-C4-18-180	502005	Oval jaws, coated, paediatric	1,8 mm	180 cm	2,0 mm	10
	BI01-C4-18-260	502006	Oval jaws, coated, balloon enteroscopy	1,8 mm	260 cm	2,0 mm	10
	BI01-C4-23-180	501900	Oval jaws, mono pull wire, coated	2,3 mm	180 cm	2,8 mm	10
	BI01-C4-23-230	501901	Oval jaws, mono pull wire, coated	2,3 mm	230 cm	2,8 mm	10
	BI01-D4-23-180	501898	Oval jaws, mono pull wire, coated, spike	2,3 mm	180 cm	2,8 mm	10
	BI01-D4-23-230	501899	Oval jaws, mono pull wire, coated, spike	2,3 mm	230 cm	2,8 mm	10
	BI01-C3-23-180	500160	Oval jaws, mono pull wire	2,3 mm	180 cm	2,8 mm	10
	BI01-C3-23-230	500161	Oval jaws, mono pull wire	2,3 mm	230 cm	2,8 mm	10
	BI01-D3-23-180	500164	Oval jaws, mono pull wire, spike	2,3 mm	180 cm	2,8 mm	10
	BI01-D3-23-230	500165	Oval jaws, mono pull wire, spike	2,3 mm	230 cm	2,8 mm	10
	BI01-G4-23-180	502162	Alligator, mono pull wire, coated	2,3 mm	180 cm	2,8 mm	10
	BI01-G4-23-230	502163	Alligator, mono pull wire, coated	2,3 mm	230 cm	2,8 mm	10
	BI01-H4-23-180	502556	Alligator, mono pull wire, coated, spike	2,3 mm	180 cm	2,8 mm	10
	BI01-H4-23-230	502557	Alligator, mono pull wire, coated, spike	2,3 mm	230 cm	2,8 mm	10

EndoBite² with mono pull wire

Ensures large sizes of sample material ideal for histological assessment because the cut surfaces are completely closed, even if the position of the endoscope is challenging. The mono pull wire prevents the branches from tilting and thus protects the working channel against wear and damage.



Biopsy

Biopsy forceps (Macro forceps)




	Product Code	Order No.	Characteristics	ø Jaws	Length	ø Working Channel	Unit
	AMHBFJ2.8*2300*	AMHBFJ2.8*2300	Alligator, mono pull wire, coated, macro forceps	2,8 mm	230 cm	3,2 mm	10
	AMHBFK2.8*2300*	AMHBFK2.8*2300	Alligator, mono pull wire, coated, macro forceps, spike	2,8 mm	230 cm	3,2 mm	10

* Only available in Europe

Hemostasis & Ligation

NEW Clipmaster® Clip applicator





	Product Code	Order No.	Characteristics	Opening Width	Length	ø Working Channel	Unit
	SPE1-X9-26-230*	502700	Coated, 360° positioning	12 mm	210 cm	2,8 mm	10

* Only available in Europe


Clipmaster® Clip applicators



	Product Code	Order No.	Characteristics	Opening Width	Length	ø Working Channel	Unit
	SPE1-X1-26-220*	502011	Clipmaster.1 1 Clip preloaded, 360° positioning	12 mm	206 cm	2,8 mm	2
	SPE1-X3-26-220*	502063	Clipmaster.3 3 Clips preloaded, 360° positioning	12 mm	206 cm	2,8 mm	2

* Only available in Europe

Balloon Enteroscopy




Compacted Wire Technology (CWT)

Compacted Wire Technology combines the benefits of a monofilament snare with those of a braided snare:

- ✓ Easier positioning with flat lesions
- ✓ Sparing tissue ablation

Polypectomy snares (CWT)



	Product Code	Order No.	Characteristics	ø Snare	Length	ø Working Channel	Unit
	POL1-B1-10-23-260-OL	502529	Convex	10 mm	260 cm	2,8 mm	10
	POL1-B1-15-23-260-OL	502530	Convex	15 mm	260 cm	2,8 mm	10
	POL1-B1-20-23-260-OL	502531	Convex	20 mm	260 cm	2,8 mm	10
	POL1-B1-30-23-260-OL	501209	Convex	30 mm	260 cm	2,8 mm	10

Polypectomy




Figure: POL1-B9-15-23-220-OL

Polypectomy snares (Monofil)



	Product Code	Order No.	Characteristics	ø Snare	Length	ø Working Channel	Unit
	POL1-B9-10-23-220-OL	502526	Oval	10 mm	220 cm	2,8 mm	10
	POL1-B9-15-23-220-OL	502050	Oval	15 mm	220 cm	2,8 mm	10
	POL1-B9-22-23-220-OL	502048	Oval	22 mm	220 cm	2,8 mm	10
	POL1-B9-30-23-220-OL	502049	Oval	30 mm	220 cm	2,8 mm	10
	POL1-B3-10-23-220-OL	502029	Oval, rotatable	10 mm	220 cm	2,8 mm	10
	POL1-B3-15-23-220-OL	501215	Oval, rotatable	15 mm	220 cm	2,8 mm	10
	POL1-B3-22-23-220-OL	501216	Oval, rotatable	22 mm	220 cm	2,8 mm	10
	POL1-B3-30-23-220-OL	501217	Oval, rotatable	30 mm	220 cm	2,8 mm	10

Balloon Enteroscopy




Conscientious and reliable cleaning

Racoon cleaning brushes ensure optimal force transmission when cleaning your endoscope. The distal protection head reliably guards the working channel. The dense brush structure makes cleaning even more efficient.

RACoon Brush Cleaning brush



	Product Code	Order No.	Characteristics	ø Brush	Length	ø Working Channel	Unit
	CLE1-C2-30-16-280	500463	2 brushes	3 mm	280 cm	2,0 - 2,8 mm	10

SHEPHERD Endoscope protector



	Product Code	Order No.	Characteristics	ø Outer	ø Inner	Length	Unit
 	SHP1-U1	502512	Four colors	37 mm	4 - 10 mm	14 cm	50

Lavage catheters



	Product Code	Order No.	Characteristics	Length	ø Working Channel	Unit
	SPE1-B5-18-120	501841	LL-connector	120 cm	2,0 mm	5
	SPE1-B5-21-120	501842	LL-connector	120 cm	2,8 mm	5



Figure: BYT1-A3-20

Bite blocks



	Product Code	Order No.	Characteristics		ø Inner	Unit
	BYT1-A3-20	500309	Standard, white, with OEKO TEX® strap		20 mm	100



FLUSHKNIFE



FlushKnife **INs**

FlushKnife **Bts**

ENHANCED USABILITY IDEAL FOR ALL PHYSICIANS

One knife covers everything from marking (coagulation), dissection to arrest of bleeding. The high versatility of this knife improves the efficiency of the procedure by reducing the exchange of instruments during the procedure.

THINNER SHEATH FOR ENHANCED SUCTION PERFORMANCE¹

EXCELLENT MANOEUVRABILITY AND STABILITY²

FLUID SUPPLY CAPABILITY WITH WATER PUMP JW-3

1. Due to tapered catheter design between distal section and handle to increase working channel lumen
2. Due to increased catheter diameter on the distal section to secure stable position in the working channel

FLUSHKNIFE BTS



FlushKnife BTS has a ball tip which touches a wider range of tissue and supports effective coagulation.

Model	Material Code	Characteristics	Identifier	Colour code	Knife length (mm)	Working length (mm)	Ø Working channel (mm)
DK2620JI	16659607	Ball tip, Slim Type	- B15 -	Green	1.5	2,000	2.8 or more
	16659619		- B20 -	Blue	2.0		
	16659621		- B25 -	Yellow	2.5		
	16659633		- B30 -	Black	3.0		
DK2623JI	16659669	Ball tip, Long slim type	- B15 -	Green	1.5	2,300	2.8 or more
	16659671		- B20 -	Blue	2.0		

FLUSHKNIFE NS



FlushKnife NS has a slim needle shaped tip.

Model	Material Code	Characteristics	Identifier	Colour code	Knife length (mm)	Working length (mm)	Ø Working channel (mm)
DK2620JI	16659554	Needle tip, Slim Type	- N10 -	White	1.0	2,000	2.8 or more
	16659566		- N15 -	Green	1.5		
	16659578		- N20 -	Blue	2.0		
	16659580		- N25 -	Yellow	2.5		
	16659592		- N30 -	Black	3.0		
DK2623JI	16659645	Needle tip, Long slim type	- N15 -	Green	1.5	2,300	2.8 or more
	16659657		- N20 -	Blue	2.0		



1. Marking (coagulation)



2. Water jet



3. Mucosal incision



4. Submucosal dissection



5. Submucosal supply



6. Coagulation

Recommended applications

Diameter	1 mm	1.5mm	2mm	2.5mm	3mm
Oesophagus	○	◎	◎	△	△
Stomach	○	◎	◎	◎	○
Colon	○	◎	◎	△	△

◎ Best indication ○ Possible Use △ Indicated in certain cases

Please note: Examples of the suitable protruding lengths are suggested by Takashi Toyonaga M.D. of Kobe University Hospital. A physician must take into consideration each particular condition of the area or lesion to be dissected when selecting a protruding knife length.

PuraStat®

Clear haemostatic control ⁽¹⁾



KEY FEATURES

- Transparent soft gel
- No preparation required
- Non-swelling ⁽²⁾
- Synthetic peptide
- Biocompatible ⁽³⁾

Mode of Action

PuraStat is a slightly viscous solution of synthetic peptides. Contact between PuraStat and blood causes the acidic peptide solution to be neutralized and exposed to ions, resulting in the formation of β -sheets that then form a 3-dimensional scaffold structure. PuraStat provides a physical barrier to stop bleeding in a variety of surgical indications.





Synthetic peptide solution

Neutral (pH) (blood, ions)

Self-assembly formation

3-D nanofiber structure

PuraStat hydrogel



Microscopic

Macroscopic

PuraStat provides a physical barrier to stop bleeding. This haemostatic agent is bioresorbable and can be used in a variety of surgical indications.

References:

- (1) PuraStat IFU-007 Rev 2
- (2) Data on file (Swell Report 2018, Eun Seok Gil 12/06/2018)
- (3) Biological Safety Report - current version

PuraStat is a class III medical device, CE marked according to European council directive 93/42/EEC on medical devices and its relatives

PuraStat®

Clear haemostatic control ⁽¹⁾

Indication for Use ⁽¹⁾

PuraStat is indicated for haemostasis in the following situations encountered during surgery, when haemostasis by ligation or standard means is insufficient or impractical:

Bleeding from small blood vessels and oozing from capillaries of the parenchyma and surrounding tissues of solid organs Examples • Myomectomy • Nephrectomy • Pancreactectomy • Partial hepatectomy • Prostatectomy • Splenectomy	Oozing from vascular anastomoses to native or artificial vessels, on the surface of blood vessels and surrounding tissues Examples • Coronary bypass • Femoral bypass • Surgery of the aorta or any peripheral arteries	Bleeding from small blood vessels and oozing from capillaries of the GI tract following surgical procedures and surrounding tissues Examples • Endoscopic mucosal resection (EMR) of GI tract • Endoscopic submucosal dissection (ESD) of GI tract • Laparoscopic resection of GI tract organs
---	--	---

PuraStat is also indicated for the **reduction of delayed bleeding** following gastrointestinal endoscopic submucosal dissection (ESD) procedures in the colon.

Please read carefully the Instructions for Use for indications and contra-indications.

Available in 1 ml (621-034), 3 ml (621-035) and 5 ml (621-036)

Advantages



Transparent Haemostat

- Maintains clear sight of bleeding point



Reapplicable & Removable

- Can be used repetitively during the same procedure
- Can be removed easily if desired



Synthetic Peptide

- Inert material
- Well accepted by the body




Easy to Use

- Single prefilled, ready-to-use syringe
- No preparation required
- Does not obstruct the applicator ⁽⁴⁾




Suitability

- Suitable for endoscopic/laparoscopic use
- Possible to place sutures, cauterize
- Applicable to narrow spaces
- Covers well uneven surfaces



B1 Droplet Reduction Mouthpiece

In response to the global pandemic, Fujifilm has identified a need for greater protection for patients and doctors alike. Fujifilm has engineered a new droplet reduction mouthpiece to reduce droplet spread. The B1 Droplet Reduction Mouthpiece was created to help reduce healthcare workers' risk of infection and to reduce patient discomfort of covering the face with a shield.



Face shield

Designed to minimize the spread of droplets from gap between patient's mouth and mouthpiece

Mouthpiece

Has two foam layers to reduce droplets from oral cavity

Accommodates different scopes

Compatible with scope diameters from 3 to 16 mm

Stable

Latex free silicone rubber band to fix the drape to patient

Designed to reduce the spread of droplets from mouthpiece.

These images illustrate a person coughing with a normal mouthpiece on and with the B1 and the dispersement of droplets from the mouth.



Normal mouthpiece



With B1

Ordering Information

Material #:	16673275
Description:	B1-DROPLET REDUCTION MOUTHPIECE W DRAPE
Material Group:	698



P1 Droplet Reduction Endoscope Cover

Fujifilm recognizes the need for greater protection of both patient and clinician during this global pandemic. Fujifilm has engineered a new droplet reduction endoscope cover to help mitigate the risk of the spread of infection. The P1 Droplet Reduction Endoscope Cover reduces the amount of droplets from device use in the working channel, protecting healthcare workers from patients' bodily fluid.



Drape Design

Covers up to the button operation portion, reducing droplets from the buttons, further reducing risk of infection.

Drape Slit Structure

Reduces discomfort during angle operation, with the same direct angle operation doctors are used to feeling during procedures.

Droplet Reduction

Main body with sponge rubber to reduce droplets when removing devices from the working channel

Compatible with most scopes

Compatible with Fujifilm G7 and G5 grips as well as competitor scopes by combining with FOV-DV7 or FV-001 forceps valves.

Designed to reduce the spread of droplets that may contain pathogens.

These images illustrate the removal of a device without the P1 Droplet Reduction Endoscope Cover and the removal of device with the P1 Droplet Reduction Endoscope Cover employed, clearly reducing the amount of droplets.



Without P1



With P1

Ordering Information

Material #:	16673316
Description:	P1-DROPLET REDUCTION ENDOSCOPE COVER
Material Group:	698

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231	Калининград (4012)72-03-81	Омск (3812)21-46-40	Сыктывкар (8212)25-95-17
Ангарск (3955)60-70-56	Калуга (4842)92-23-67	Орел (4862)44-53-42	Тамбов (4752)50-40-97
Архангельск (8182)63-90-72	Кемерово (3842)65-04-62	Оренбург (3532)37-68-04	Тверь (4822)63-31-35
Астрахань (8512)99-46-04	Киров (8332)68-02-04	Пенза (8412)22-31-16	Тольятти (8482)63-91-07
Барнаул (3852)73-04-60	Коломна (4966)23-41-49	Петрозаводск (8142)55-98-37	Томск (3822)98-41-53
Белгород (4722)40-23-64	Кострома (4942)77-07-48	Псков (8112)59-10-37	Тула (4872)33-79-87
Благовещенск (4162)22-76-07	Краснодар (861)203-40-90	Пермь (342)205-81-47	Тюмень (3452)66-21-18
Брянск (4832)59-03-52	Красноярск (391)204-63-61	Ростов-на-Дону (863)308-18-15	Ульяновск (8422)24-23-59
Владивосток (423)249-28-31	Курск (4712)77-13-04	Рязань (4912)46-61-64	Улан-Удэ (3012)59-97-51
Владикавказ (8672)28-90-48	Курган (3522)50-90-47	Самара (846)206-03-16	Уфа (347)229-48-12
Владимир (4922)49-43-18	Липецк (4742)52-20-81	Саранск (8342)22-96-24	Хабаровск (4212)92-98-04
Волгоград (844)278-03-48	Магнитогорск (3519)55-03-13	Санкт-Петербург (812)309-46-40	Чебоксары (8352)28-53-07
Вологда (8172)26-41-59	Москва (495)268-04-70	Саратов (845)249-38-78	Челябинск (351)202-03-61
Воронеж (473)204-51-73	Мурманск (8152)59-64-93	Севастополь (8692)22-31-93	Череповец (8202)49-02-64
Екатеринбург (343)384-55-89	Набережные Челны (8552)20-53-41	Симферополь (3652)67-13-56	Чита (3022)38-34-83
Иваново (4932)77-34-06	Нижний Новгород (831)429-08-12	Смоленск (4812)29-41-54	Якутск (4112)23-90-97
Ижевск (3412)26-03-58	Новокузнецк (3843)20-46-81	Сочи (862)225-72-31	Ярославль (4852)69-52-93
Иркутск (395)279-98-46	Ноябрьск (3496)41-32-12	Ставрополь (8652)20-65-13	
Казань (843)206-01-48	Новосибирск (383)227-86-73	Сургут (3462)77-98-35	
Россия +7(495)268-04-70	Киргизия +996(312)-96-26-47	Казахстан +7(7172)727-132	